
FastGene™ RNA Premium Kit

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FastGene™ RNA Premium Kit

For purification of total RNA (gDNA removed) from animal cells/tissues

FG-81006 (6preps), FG-81050 (50preps), FG-81250 (250preps)

Step	Standard protocol	Large input protocol
Sample quantity	< 5×10 ⁶ cultured animal cells < 10 mg animal tissues	< 1×10 ⁷ cultured animal cells < 20 mg animal tissues
Resuspension/homogenization by cell lysis	 350 μL buffer RL ^{※1} Vortex vigorously 	 600 μL buffer RL ^{※1} Vortex vigorously 
Filtration of cellular debris	 Transfer lysate into a FastGene™ RNA filter column Centrifuge at ≥ 10,000 x g (RT : 20 ~ 25°C) 1 min Discard filter column, Harvest flow-through	
Optimize RNA binding conditions	 Add 350 μL 70 % ethanol Mix thoroughly by pipetting	 Add 600 μL 70 % ethanol Mix thoroughly by pipetting
RNA binding	 Load mix (up to 700 μL) onto FastGene™ RNA binding column Centrifuge at ≥ 10,000 x g (RT : 20 ~ 25°C) 1 min Discard flow-through, Re-insert binding column in collection tube (2.0 mL)	Repeat the procedure for larger volume.
Protein elimination	 Add 600 μL of buffer RW1 Centrifuge at ≥ 10,000 x g (RT : 20 ~ 25°C) 30 s Transfer binding column to new collection tube (2.0 mL)	
Desalination	 Add 700 μL of buffer RW2 ^{※1} Centrifuge at ≥ 10,000 x g (RT : 20 ~ 25°C) 30 s Transfer binding column to new collection tube (2.0 mL)	
Removal of RW2	 Centrifuge at full speed (RT : 20 ~ 25°C) 1 min Transfer binding column to new collection tube (1.5 mL)	
Elution of RNA	 Add 50 μL of buffer RE to membrane center Centrifuge at ≥ 10,000 x g (RT : 20 ~ 25°C) 1 min Discard binding column, Harvest eluted solution 	
Optimize DNase I conditions	 Add 5 μL 10x DNase I reaction buffer	
DNA digestion	 Add 1 μL of DNase I ^{※1} to the mixture Mix thoroughly by pipetting Incubate for 10 min (at room temp. 20-25°C)	
RNA rebinding optimization	 Add 250 μL of buffer RBD ^{※1} to the mixture Mix thoroughly by pipetting	
RNA binding	 Transfer mixture into FastGene™ RNA mini-elute column Centrifuge at ≥ 10,000 x g (RT : 20 ~ 25°C) 1 min Discard flow-through, Re-insert binding column in collection tube (2.0 mL)	
Desalination/ Elimination of digested DNA	 Add 700 μL buffer RW2 ^{※1} Centrifuge at ≥ 10,000 x g (RT : 20 ~ 25°C) 30 s Transfer mini-elute column in new collection tube (2 mL)	
Removal of RW2	 Centrifuge at full speed (RT : 20 ~ 25°C) 1 min Transfer mini-elute column in new collection tube (1.5 mL)	
Elution of RNA	 Add appropriate volume ^{※2} of buffer RE to the membrane center Centrifuge at ≥ 10,000 x g (RT : 20 ~ 25°C) 1 min Discard mini-elute column, Harvest eluted solution	

※1 : need preparation before use. ※2 : Standard protocol 20 μL (10 ~ 50 μL), Large input protocol 50 μL (20 ~ 50 μL)

 Safety Stopping Point. Sample after homogenize step or 1st elution step can be stored at -70°C.

NIPPON GENETICS EUROPE GmbH

(English web page) www.nippongenetics.eu (German web page) www.nippongenetics.de

2017 NOV/Rev.01.0