# KAPATaq™ EXtra DNA ポリメラーゼ

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Cat.No	Kit Size
<b>KK3009</b>	<b>250units</b>
<b>KK3008</b>	<b>50units</b>

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Tel. 075 (257) 5421 Fax. 075 (257) 5422

## Quick Reference Guide

Product code	Kit size
<b>KK3008</b>	<b>50 units</b>
<b>KK3009</b>	<b>250 units</b>

### Kit Components

- 5 U/μl KAPATaq EXtra DNA Polymerase
- 5x KAPATaq EXtra Buffer
- 25 mM MgCl<sub>2</sub> solution
- KAPA dNTP Mix (10mM each)

### Storage

Store all components at -20 °C.

## Production Description

The KAPATaq EXtra PCR system is a blend of *Taq* DNA polymerase and a modified archaeal (Type B) DNA polymerase possessing proofreading capability. This two-enzyme system is designed specifically to support robust, long range, and sensitive PCR. The KAPATaq EXtra system polymerizes DNA from a primer annealed to a DNA template in the presence of deoxyribonucleotide triphosphates. Both enzymes possess 5'-3' polymerase activity, but only *Taq* possesses 5'-3' exonuclease activity and only the Type B polymerase possesses 3'-5' exonuclease (proofreading) activity.

KAPATaq EXtra displays higher fidelity than *Taq* polymerase.

KAPAHifi™ DNA Polymerase is recommended for shorter targets where higher fidelity is required.

## Applications

The KAPATaq EXtra PCR system is suited for:

- PCR amplification of long targets and/or PCR using low concentrations of template DNA
- Standard short- and mid-range PCR amplification
- Production of PCR products to be used for ligation into 3'-T-overhang cloning vectors

## KAPATaq EXtra PCR Protocols

### Standard PCR Protocol for short targets (up to 8 kb) and/or high concentrations of template DNA

Typically, the PCR reaction setup might consist of:

	Final concentration	50μl rxn
PCR grade water up to 50μl		As required
5x KAPATaq EXtra Buffer (without Mg <sup>2+</sup> )	1x	10μl
MgCl <sub>2</sub> (25mM)	{ 1.75mM 2.5mM	3.5μl 5.0μl *
dNTPs (10mM each dNTP)	0.3mM	1.5μl
Fwd primer (10μM)	0.5μM	2.5μl
Rev primer (10μM)	0.5μM	2.5μl
Template	As required	As required
KAPATaq EXtra DNA Polymerase (5U/μl)	1.25U/50μl	0.25μl
<b>Total</b>		<b>50μl</b>

\*Please try two MgCl<sub>2</sub> concentrations.

### PCR cycling conditions might consist of:

Initial Denaturation:	94°C	2min	} 25 cycles**
Denaturation:	94°C	15 - 25sec*	
Annealing:	T <sub>m</sub> -5°C	15 sec	
Extension:	72°C**	1min per 1kb	
Final Extension:	72°C	1min per 1kb	

\* 15 sec for slow-block instruments, 25 sec for fast-block instruments.

\*\*For longer targets (up to to 12 kb) and/or lower concentrations of template DNA lower the extension temperature to 68°C and perform 35 cycles.

## License

The purchase of this product conveys to the purchaser only the limited, non-transferable right to use the purchased quantity of the product for the purchaser's own research by the purchaser only under the following U.S. patent claims and foreign counterpart patent claims: U.S. Patent No. 5,436,149 (claims 6-16). No rights are granted to the purchaser to sell, modify for resale or otherwise transfer this product, either alone or as a component of another product, to any third party. Takara Bio reserves all other rights, and this product may not be used in any manner other than as provided herein. For information on obtaining a license to use this product for purposes other than research, please contact Takara Bio Inc., Seta 3-4-1, Otsu, Shiga 520-2193, Japan (Fax +81-77-543-9254).