Taiwan
Advanced
Nanotech

TANBead Nucleic Acid Extraction
System (non-sterile)

039.M0101.X01 Rev V1.0_2021.01

Maelstrom 8 Autostage User Manual

Table of Contents

1.	Introduction	4
2.	Instrument Overview	9
3.	Before Usage	12
4.	Attention of Installing Spin Tips	13
5.	Home Page	14
6.	Guide	17
7.	Manual	18
8.	Edit	20
9.	Setting	28
10.	Error Message	31
11.	Cleaning & Maintenance	33
12.	Disposal	33
13.	Patent Info	34
14.	Manufacturer	34

About this manual

The label on the instrument, the User Manual, and other packaging material may contain the following symbols.

Catalog number

Serial number

Specification of fuse

Manufacturer

The date of manufacture

This product fulfills the requirements of the European Directive

In Vitro Diagnostic Medical Device

Consult Instructions for Use

Authorized Representative in the European Community

Caution, consult accompanying documents

Hot surface, contact with skin may cause burns

Watch your fingers and your hands

1. Introduction

Magnetic rotary mixer is a magnetic bead handling device with high speed stirring function, which can spin up to 3000 rpm. With eight magnetic rods, intuitive interface and simple operation, it can accomplish any application of magnetics beads. It works alone with Autostage to form an automated solution, called Maelstrom 8 Autostage.

Operation Principle

Maelstrom 8 Autostage uses magnetic rods to collect and transfer magnetic beads and spin tips to mix suspension. Rotating the spin tips will mix suspension very well. Extended magnetic rods can collect magnetic beads efficiently. Purified nucleic acids can be obtained after cell lysis, nucleic acids absorption, wash and elution.

Intended Use

Maelstrom 8 Autostage can be used in combination with TANBead Nucleic Acid Extraction Kit to extract DNA/RNA from specimen.

This medical device is used by physicians or medical technologists.

1. Introduction

Environmental Requirements

- Choose a location with good air circulation
- Place this instrument on a table that can bear the weight 15 kg.
- Do not use this instrument in a location where there is large temperature variability or high humidity.
- Operation condition :
Temperature : 4-40 °C
Humidity : 40-80 %
- Storage and transportation condition :
Temperature : 8-40 °C
Humidity : 0-85 %
- Maximum operating altitude : 2000 m

1. Introduction

Safety instructions and guidelines

- This device can be used with potentially biohazardous materials. Use appropriate personal protective equipment (gloves, safety goggles, lab coat, etc.) for handling and disposing of biohazardous materials.
- Under normal condition, sound pressure level from Maelstrom 8 Autostage does not exceed 65 dB and cause a hazard. Please contact technical support in case of higher sound pressure level.
- This device can be hazardous due to the use of chemical and biohazardous substances.
- Users should adhere to their institutional guidelines for the handling and disposal of all infectious substances used with this device.
- It is important to clean the device after every use. If samples or reagents have been spilled, it is important to clean the instrument immediately to avoid damage or contamination of samples.
- This device is to use with the compatible spin tips. Using incompatible spin tips may cause poor extraction performance.
- Read this user manual in its entirety prior to operating the device. Failure to read, understand, and follow the instructions in the manual may result in damage to the device, injury to laboratory and operating personnel or poor performance.

1. Introduction

ENGLISH

Safety Requirements

- The medical device has passed the tests and conformed to the standards of IEC 61010-1:2010, "Safety requirements for electrical equipment for measurement, control, and laboratory use - Part 1: General requirements".
- The medical device has passed the tests and conformed to the standards of IEC 61010-2-101:2015, "Safety requirements for electrical equipment for measurement, control and laboratory use - Part 2-101: Particular requirements for in vitro diagnostic (IVD) medical equipment".

EMC Requirements

- The medical device has passed the tests and conformed to the standards of IEC 61326-1:2012 & IEC 61326-2-6:2012, "Electrical equipment for measurement, control and laboratory use - EMC requirements - Part 2-6: Particular requirements - In vitro diagnostic (IVD) medical equipment".
- The medical device has passed the tests and conformed to following standards:
 - EN 61326-1:2013
 - EN 61326-2-6:2013

1. Introduction

What's in the Box

Please note that following accessories accompanying with instrument may vary from region to region.

Magnetic rotary mixer x1

Spin tips x8

User manual x1

Autostage x1

Power cord & adapter x1

Transmission cable x1

Power cord plug may vary from region to region

2. Instrument Overview

ENGLISH

2. Instrument Overview

ITEM	SPECIFICATION
Model Name	Magnetic rotary mixer
Weight	600 g
Dimensions	112(W)x327(H)x63(D)mm
Power	5 Vac, 2 A
Battery	3.7 Vac , 2850mAh
Throughput	8 samples per run
Pro. volume	50 μ l ~ 1500 μ l
Spin Speed	up to 3000 rpm
Magnetic rod	>3000 gauss
Display	2.4" LCD, 240 x 320 pixels



2. Instrument Overview

ENGLISH

ITEM	SPECIFICATION
Model	Autostage
Weight	9.6 kg
Dimensions	*357(W)x285(H)x190(D) mm
Power	19 Vac, 120 W
Heater	1 heating block

* Without cover open

3. Before Usage

Each magnetic rotary mixer is set into shipping mode before it's delivered, there's a need to unlock the shipping mode before first use. Once it is plugged in, the shipping mode will be unlocked.

Unlock Shipping Mode & Set up Maelstrom 8 Autostage

- Connect the magnetic rotary mixer with Autostage via transmission cable
- Plug the Autostage to the power (19V adapter required)
- Power on magnetic rotary mixer to unlock shipping mode, and followed by Autostage
- Go to Setting\System Param Setting to enable & apply "Stand Connection"

Transmission Cable

- Side I: Micro USB and Mini USB
 - Micro USB to charging port of magnetic rotary mixer
 - Mini USB to Mini USB port of magnetic rotary mixer
- Side II: Mini USB and Type A USB
 - Mini USB to Autostage communication port
 - Type A USB to Autostage Type A USB port

Magnetic Rotary Mixer



- When the rotary mixer is being charged, the charging indicator will turn orange.
- **Battery power will decrease slowly even the rotary mixer is not being used.**
- DO NOT disassemble the magnetic rotary mixer
- Pressing MENU to power on or turn off the device

4. Attention of Installing Spin Tips


When installing spin tips on magnetic rotary mixer, please pay attention to following requirement.

- When mount the tips on magnetic rotary mixer, make sure the leading edge of spin shaft is as close as possible to the bottom of sleeve. Too much interspace will lead to abnormal operation.

5. Home Page

From Main page, user can choose Guide mode or Manual mode to operate the rotary mixer. To edit programs, please enter Edit page. To set up system parameter, such as beep sound, auto-off, time and date, please go to Setting page. Further system information can be found in Setting page.

5. Home Page


Functional Keys – MENU button

- 1) Press MENU button to power on the magnetic rotary mixer
- 2) Press MENU button for 3 seconds to power off the magnetic rotary mixer
- 3) During the operation, press MENU button once to interrupt the experiment and then follow the prompt message to press second time to cease the operation.
- 4) User can return to main page without saving changes by pressing MENU button twice.

Functional Keys – ROD button

- 1) In Manual mode, press ROD button to extend or withdraw the magnetic rods.
- 2) In Edit page, press ROD button to switch the input method

5. Home Page


Functional Keys – ENTER button

- 1) In Guide mode, press ENTER to start the protocol selected.
- 2) During the extraction process, press ENTER once to pause the operation. Then, press ENTER second time to resume the operation.
- 3) In Manual mode, press ENTER to start the mixing.
- 4) In Edit page, in each field, press ENTER to start to modify the parameter value of field.
- 5) To save the parameter value, use click wheel to move the cursor to SAVE and press ENTER to save.

Functional Keys – Click Wheel


- 1) Click wheel is mainly used to move cursor to different fields of Edit page.
- 2) In Manual mode, press up/down of click wheel to speed up/down the spin speed.
- 3) In Edit page, press up/down of click wheel to adjust the parameter value set in the field.

6. Guide

Pre-installed Programs

On each magnetic rotary mixer, there are several programs pre-installed before it is shipped from the manufacturing plant. Use the click wheel to choose desired program for the experiment.

Know the Icons

In the course of experiment, the display screen will automatically switch icons to represent current state of experiment. Knowing those icons is important.

Mixing

Beads
Collecting

Vapor

Finish of
Protocol

7. Manual

For the operation without using Autostage

User can do the nucleic acid extraction without using Autostage. In this case, manual mode shall be applied.

How to do the extraction manually

For manually extraction, there is no need to use Autostage.

Mixing – Manual Operation

- 1) Mount the tips to magnetic rotary mixer (tips combo pack required)
- 2) Insert it to the 96 deep well plate and press ENTER button to start mixing. User can press up/down of click wheel to accelerate or decelerate the spin speed. In addition, user can press left/right of click wheel to alter the spin direction.
- 3) Press the ENTER button to cease mixing

7. Manual

For the operation without using Autostage


User can do the nucleic acid extraction without using Autostage. In this case, manual mode shall be applied.

How to do the extraction manually

For manually extraction, there is no need to use Autostage.

Beads Collecting – Manual Operation

- 1) Raise magnetic rotary mixer off the plate and press ROD button to extend magnetic rods
- 2) Slowly move down to plate well and beads will be collected on the leading edge of spin tips
- 3) Raise the mixer off the plate and slowly move to targeted plate wells
- 4) Submerge the spin tips with collected beads in the solution buffer and then press ROD button to withdraw magnetic rods. After withdraw rods, beads will be released in the solution buffer. Then, repeat mixing for binding.


8. Edit

Editing Protocol


There are some pre-installed programs on magnetic rotary mixer. User can freely edit them on either rotary mixer or by importing a program file.

How to import a program file to magnetic rotary mixer

1. Link the rotary mixer with the PC by the transmission cable.
 - » Type A USB to the PC
 - » Micro USB to the charging port of magnetic rotary mixer

2. Power on the rotary mixer
3. Go to Windows Explorer and navigate Computer, then find "USB Drive" where you will see all pre-installed files.

4. Copy and paste targeted file (*.csv) to the USB Drive

 The *.csv file name MUST be a 2-digit number (ex. 01~50)

8. Edit

Protocol File

Please note that the protocol is a *.csv file. User can copy pre-installed protocol file and paste to a PC for further editing use.

	A	B	C	D	E	F	G	H	I
1	20190910-001								
2		Well:1	Well:2	Well:3	Well:4	Well:5	Well:6	Well:7	Well:8
3	Well capacity	1000	1000	900	800	900	50	900	9
4	Step	Well	Action (C:Colle	RPM	Time(Second)	CW/CCW(Secon	Temperature	Temperature_Control	
5	1	2	Mixing	2500	10	0	35	NO	
6	2	2	Collection	0	10	0	35	NO	
7	3	2	Pause	0	10	0	35	NO	
8	4	1	Mixing	2000	10	0	60	No	
9	5	1	Collection	0	10	0	35	NO	
10	6	3	Mixing	2500	10	0	50	NO	
11	7	3	Collection	0	10	0	50	NO	
12	8	3	Pause	0	10	0	45	NO	
13	9	4	Mixing	2500	10	0	20	NO	
14	10	4	Collection	0	10	0	20	NO	
15	11	5	Mixing	2000	15	0	20	NO	
16	12	5	Collection	0	15	0	20	NO	
17	13	5	Vapor	0	10	0	20	NO	
18	14	6	Mixing	1800	15	0	50	No	
19	15	6	Collection	0	10	0	20	NO	
20	16	5	Mixing	2000	10	0	20	NO	
21									

Know the Protocol

Program Name

- Stored in A1 field of spreadsheet
- Maximum 14 characters

Well capacity (=Volume in the well)

- Input value from 50µl to 1500µl
- Well 1~12 (B3 to M3 field of spreadsheet) all need a value

Step

- Maximum 80 steps in a program

Well

- Choose a well to set specific action
- Maximum 12 wells capable of being proceeded

8. Edit

Know the Protocol

Action

- Give specific requirement to a well for extraction, such as mixing, collection, pause, or vapor

RPM

- Speed of spin
- Input value from 500 to 3000

Time

- Operation time of each action
- Input value on the basis of second
- Maximum 3599 seconds

CW/CCW (Clockwise/Counterclockwise)

- Frequency of changing spin direction, ex. 20 means spin direction will change every 20 seconds
- Input value on the basis of second, from 0~999

Temperature

- Temperature requirement of each well for action
- Input value on the basis of centigrade, from 20-110C
- Increase/decrease in increments of 5 degrees

Temperature Control

- Input value YES or NO
- YES means heating block needs to reach targeted temp before mixing, collection, or vapor
- NO means temperature has nothing to do with action of well

8. Edit

How to Edit Program on Magnetic Rotary Mixer


Editing protocol parameter needs the use of click wheel.

Procedure of Edit


1. Go to Edit page
2. Select a program to edit
3. Choose "Edit" to apply
4. Use left/right key to move to desired field or use up/down key to switch page
5. Select a field to change input value and then press ENTER button to access
6. Press up/down to change value in the field and then press ENTER to confirm
7. Once the adjustment has been done, be sure to SAVE the adjustment


Go to Edit page


Select a program


Choose Edit function

Left/Right key to move to desired field

ENTER to access for adjustment


Up/Down key to adjust and ENTER to confirm


8. Edit

How to Edit Program on Magnetic Rotary Mixer

Editing protocol parameter needs the use of click wheel.

After value is adjusted, press ENTER to confirm and the field background will turn to Yellow.

When the adjustment is complete, be sure to move to SAVE icon and press ENTER to memorize all the updated value.

Change editing page

8. Edit

About Protocol Setting

Knowing the meaning of each field in Edit page is important.



14:11 100% [Battery Icon]

Cap. & Temp. Edit

01	0700	02	0700	03	0800
04	0800	05	0800	06	0150
07	0700	08	0800	09	0800
10	0800	11	0800	12	0150

Temperature: 000 Enable:

SAVE << >>

The temperature requirement here is a rapid & apply-to-all setting. Once the value is given and Enable is checked, the value will be applied to all plate wells.

8. Edit

About Protocol Setting

Knowing the meaning of each field in Edit page is important.

14:11 100% [Battery Icon]

Edit

03 : 611-1

Step : 01 Mode : Mixing

Well : 0004 Time : 01:00

Stir : 000 RPM : 3000

Temperature 050 Enable V

SAVE << >>

03: 611-1

Program name

Number of program

Step: 01 Mode: Mixing

Order of step

Action: Mixing, Collection, Pause, or Vapor

Well: 0004 Time: 01:00

Number of well

Time of executing action
Value from 00:01-59:59

Stir: 000 RPM: 3000

Spin speed, value from 500-3000


Frequency of changing spin direction
Value from 0-999 seconds

8. Edit

About Protocol Setting

Knowing the meaning of each field in Edit page is important.


14:11 100% [|||||]
Edit
03 : 611-1
Step : 01 Mode : Mixing
Well : 0004 Time : 01:00
Stir : 000 RPM : 3000
Temperature 050 Enable
SAVE << >>

The temperature requirement here is a setting applied to defined plate well only.

9. Setting

In Setting page, the user can find system info and some basic system parameters.

There are some basic settings, such as Time, Auto-off, Beep & LCD, and System Parameter, which the user can freely change.

SW Version:

Magnetic rotary mixer software version

HW Version:

Magnetic rotary mixer hardware version

S/N:

Serial number of magnetic rotary mixer

9. Setting

In Setting page, the user can find system info and some basic system parameters.

Time Setting

Use click wheel to set up correct time of device. Once the setup is done, move to Apply to save changes.

Auto-off Setting

Time setting of Auto-sleep and Auto-shutdown is for handler only.

When the rotary mixer is connected to Autostage and Autostage power is on, the Auto-shutdown setting will be invalid.

Shipping Mode

Each magnetic rotary mixer is set into shipping mode before it's delivered, once it is plugged in, shipping mode will be unlocked.

9. Setting

In Setting page, the user can find system info and some basic system parameters.

Beep & LCD Setting

User can turn on/off the button beep at will.

LCD Backlight: from 1 to 10

System Param Setting

It is necessary to Enable stand connection when a magnetic rotary mixer is in combination with Autostage as an automated system.

If stand connection is set in Disable mode, even the rotary mixer is connected with Autostage thru transmission cable, the nucleic acid extraction will not be made automatically.

Rod Offset

Rod offset value is set by the manufacturer and it significantly influences the performance of nucleic acid extraction. We DO NOT suggest user changing this value and we and our distributors are not responsible for any consequences coming from changed value.

10. Error Message

Understanding the meaning of error message may help you solve the problem in use.

No.	Description	Error Message	LED	Alarm
1	<p>System image file is not found.</p> <p><u>How to deal</u> Try to reboot. If still failed, please contact our technician.</p>		N/A	NO
2	<p>Battery is too low</p> <p><u>How to deal</u> Please plug in rotary mixer for charging</p>		Red	YES
3	<p>Battery voltage is less than 2V</p> <p><u>How to deal</u> Please contact our technician for further help</p>		Yellow	NO
4	<p>Battery temp is too high (over 45C)</p> <p><u>How to deal</u> Please stop using device for 10 mins and then try again.</p>		N/A	NO

10. Error Message

Understanding the meaning of error message may help you solve the problem in use.

No.	Description	Error Message	LED	Alarm
5	<p>Detecting motor module failed</p> <p><u>How to deal</u> Try to reboot. If still failed, please contact our technician.</p>	 <p>The screenshot shows the TAN Bead logo at the top, followed by the text "Warning: Cannot find Motor Module." in red at the bottom.</p>	N/A	NO
6	<p>DC motor spin speed incorrect</p> <p><u>How to deal</u> Try to reboot. If still failed, please contact our technician.</p>	 <p>The screenshot shows a yellow warning bar at the top with the text "Warning". Below it is a yellow warning icon and the text "Motor block error!".</p>	N/A	NO
7	<p>Stepping motor over current</p> <p><u>How to deal</u> Please contact our technician for further help</p>	 <p>The screenshot shows a yellow warning bar at the top with the text "Warning". Below it is a yellow warning icon and the text "Step Motor OC error!".</p>	N/A	NO
8	<p>Z-axis failed in initialization or sensor no response</p> <p><u>How to deal</u> Try to reboot. If still failed, please contact our technician.</p>	 <p>The screenshot shows the TAN Bead logo at the top, followed by the text "Warning: Limit sensor no response." in red at the bottom.</p>	N/A	NO

11. Cleaning & Maintenance

Periodically cleaning and maintenance is highly advised.

- It is important to clean both rotary mixer and Autostage after every use. If samples or reagents have been spilled, it is necessary to clean immediately to avoid damage or contamination of samples.
- Wear gloves and appropriate personal protective equipment. If the device is used with biohazardous materials, dispose of any cleaning materials used in accordance with your institutional guidelines.
- The device may go through a run with the magnetic rods unprotected. If this happens, the magnetic rod needs to be cleaned immediately.
- To clean the magnetic rods, wipe with a soft cloth dampened with pure water. Do not use alcohol solvent.
- If the magnetic rods cannot be cleaned, please contact us (Tel: +886-3-3167568) for technical assistance.

12. Disposal

In case the product is to be disposed of, the relevant legal regulations are to be observed.

As disposal regulations varies from country to country, please contact your supplier if necessary. Please must follow your institutional and country-specific requirements to handle the disposal of instrument and accessories. Device must be decontaminated prior to disposal.

13. Patent Info

Patents owned by Taiwan Advanced Nanotech Inc.

Patent List	
USA	US09616398B2
EU	EP2937136
Canada	CA2862946
Japan	JP6151735B2
Korea	10-1696517
China	CN104971638B
Taiwan	TWI526245B
WIPO	WO2016127292

14. Manufacturer

- Manufacturer : Taiwan Advanced Nanotech Inc.
- Manufacturer Add : No. 2, Aly. 12, Ln. 81, Longshou St., Taoyuan District., Taoyuan City 330, Taiwan
- Manufacturer Tel : +886-3-3607555

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