

1 and 7 before adding the liquid, and discard 40μL Elution buffer from column 5 and 11.)

5) Place 96-deep-well plate to the instrument, then plug in 8-strip tips and run the program.

Step	Name	Well	Volume (μL)	Waiting Time (min : ss)	Speed	Mixing Time (min : ss)	Adsorption	Magnet Time (min : ss)
1	Lysis	1	710	-- : --	fast	10 : 00	N	-- : --
2	Lysis	1	710	00 : 05	fast	00 : 00	N	-- : --

Lysis temperature:65°C, lysis stop heating step: 2th

6) At the end of step 5, take the 96-deep-well pre-packed plate out. Add 300μL isopropanol into column 1 and 7, and put the 96-deep-well plate back into NPA-32P nucleic acid purification system. The automatic extraction experiment was carried out according to the following procedure:

Step	Name	Well	Volume (μL)	Waiting Time (min : ss)	Speed	Mixing Time (min : ss)	Adsorption	Magnet Time (min : ss)
1	Beads	6	200	-- : --	slow	00 : 30	N	00 : 30
2	Bind	1	1000	-- : --	fast	08 : 00	S	02 : 00
3	Wash 1	2	800	-- : --	fast	03 : 00	S	01 : 00
4	Wash 2	3	600	-- : --	fast	03 : 00	S	01 : 00
5	Wash 3	4	600	-- : --	fast	03 : 00	S	01 : 00
6	Elution	5	100	05 : 00	fast	10 : 00	C	01 : 40
7	Discard	6	200	-- : --	slow	00 : 30	N	-- : --

Elution temperature: 60°C, Elution start heating step: 6th

7) After the automation process is completed, the purified nucleic acid is in column 5; if not used immediately, transfer the elution to a new 0.5 mL nuclease-free centrifuge tube and save it at -80 °C.

【Analysis Nucleic Acid】

Get some DNA, diluted in a advisable factor with elution buffer.

Survey the OD260, OD280 and OD320.

Concentration (μg/mL) = 50×OD260×dilution fact

1.7≤OD260-320/ OD280-320≤2.1

【Basic Information】

Hangzhou Bioer Technology Co.,Ltd

Address: No.1192 Bin'An Rd, Binjiang District, Hangzhou, Zhejiang Province, China

Tel: 0571-87774567 Fax: 0571-87774553

Web: www.bioer.com.cn

Zip Code: 310053

Aftersales Service Provider: Hangzhou Bioer Technology Co.,Ltd

MagaBio plus General Genomic DNA Purification Kit II

【Introduction】

MagaBio plus General Genomic DNA Purification Kit II is designed for rapid and reliable isolation of high-quality genomic DNA from blood, tissue, saliva, buccal swabs, cultured cells, and dry blood spots. MagaBio Particles, magnetic beads, provide a quick magnetic response time reducing overall processing time. This system combines the reversible nucleic acid-binding properties of MagaBio paramagnetic particles with the time-proven efficiency of Bioer automatic nucleic acid purification system to provide a fast and convenient method to isolated DNA from a variety of samples. Utilizing paramagnetic particles provides high-quality DNA that is suitable for direct use in most downstream applications, such as amplification and enzymatic reactions, next-generation sequencing.

【Kit Components】

Cat#	BSC74S1E	BSC74S1B	BSC74M1B
Components	32Tests	50Tests	100Tests
TES Buffer	16 mL	25 mL	50 mL
PK Solution	320 μL	500 μL	1 mL
RNase A Solution	32 μL	50 μL	100 μL
Lysis Buffer		20 mL	40 mL
Wash Buffer I	96 well pre-packed plate 2 pieces	16 mL (Add 24 mL absolute ethanol before use)	32 mL (Add 48 mL absolute ethanol before use)
Wash Buffer II		18 mL (Add 42 mL absolute ethanol before use)	18 mL×2 (Add 42 mL absolute ethanol before use)
Elution Buffer		10 mL	20 mL
Maganetic beads		1.5 mL	3 mL
Handbook		/	/

【Storage】

- The kit can be transported at room temperature.
- All reagents are stored at 2-8°C.
- All reagents, when stored properly, are stable for 12 months from the time of delivery.

【Principle and Advantage】

Nucleic acid is released under the action of Proteinase K Solution and Lysis Buffer. Magnetic beads can combine with DNA exclusively and specifically after adding isopropanol. Magnetic beads combined with DNA are captured by magnetic rod; Multiple Wash Buffers will remove contaminants. Purified DNA is released into the Elution Buffer by Magnetic beads.

MagaBio technic has great advantages:

- Mini sample, high purification;
- Simple and streamline separation procedure, used for auto-platform;
- First elution can acquire 85% or more DNA;
- No high salt solution. No inhibitor;
- No spin column.

【Materials and Equipment to be prepared by user】

1. Magnetic Rack or Bioer NPA-32P purification instrument;
2. Water bath or Dry bath;
3. Vortex mixer;
4. Absolute alcohol (For BSC74S1B and BSC74M1B).

【Reagent preparing】

1. For BSC74S1B and BSC74M1B

Add 400 μL Lysis Buffer to the column 1 and 7 of 96-deep-well plate. Add 800μL Wash Buffer I to the column 2 and 8. Add 600μL Wash Buffer II to the column 3, 4, 9 and 10. Add 100μL Elution Buffer to the column 5 and 11. Add 170μL pure water and 30μL Maganetic beads to the column 6 and 12.

2. For BSC74S1E

Mix the pre-packed plate by inverting, and centrifuge the 96-well plate in plate centrifuge for a few seconds. Tear off aluminum foil film of 96-deep-well pre-packed plate. Identify the direction of the plate (magnetic beads is in column 6th&12th).

【Protocol】

1. Blood or saliva processing

- 1) Add 200-250μL blood or saliva and 10μL PK Solution to the column 1and 7 of 96-deep-well plate. Avoid cross-contamination.
- 2) Place 96-deep-well plate to the instrument, then plug in the 8-strip tips and run the program.

Step	Name	Well	Volume (μL)	Waiting Time (min : ss)	Speed	Mixing Time (min : ss)	Adsorption	Magnet Time (min : ss)
1	Lysis	1	710	-- : --	fast	10 : 00	N	-- : --
2	Lysis	1	710	00 : 05	fast	00 : 00	N	-- : --

Lysis temperature:65°C, lysis stop heating step: 2th

- 3) At the end of step 2, take the 96-deep-well pre-packed plate out. Add 300μL isopropanol into column 1 and 7, and put the 96-deep-well plate back into NPA-32P nucleic acid purification system. The automatic extraction experiment was carried out according to the following procedure:

Step	Name	Well	Volume (μL)	Waiting Time (min : ss)	Speed	Mixing Time (min : ss)	Adsorption	Magnet Time (min : ss)
1	Beads	6	200	-- : --	slow	00 : 30	N	00 : 30
2	Bind	1	1000	-- : --	fast	08 : 00	S	02 : 00
3	Wash 1	2	800	-- : --	fast	03 : 00	S	01 : 00
4	Wash 2	3	600	-- : --	fast	03 : 00	S	01 : 00
5	Wash 3	4	600	-- : --	fast	03 : 00	S	01 : 00
6	Elution	5	100	05 : 00	fast	10 : 00	C	01 : 40
7	Discard	6	200	-- : --	slow	00 : 30	N	-- : --

Elution temperature: 60°C, Elution start heating step: 6th

2. Tissue processing

- 1) Prepare all reagents and samples under room temperature. Grind the tissue into powder under the liquid nitrogen
- 2) Add 30mg sample powder to the 1.5mL DNase-free tube. Add 10μL PK Solution and 300μL TES Buffer into the tube and mix by pulse-vortexing for 15 seconds.
- 3) Incubate at 56°C for 30 minutes (The lysis step can be more longer for 1-3h,or overnight for ear and tail) .
- 4) Remove the tube from 56°C.Centrifuge for 5min at 12000g. Transfer 300μL of the supernatant to the column 1 and 7 of 96-deep-well plate.

Place 96-deep-well plate to the instrument, then plug in 8-strip tips and run the program.

Step	Name	Well	Volume (μL)	Waiting Time (min : ss)	Speed	Mixing Time (min : ss)	Adsorption	Magnet Time (min : ss)
1	Lysis	1	710	-- : --	fast	10 : 00	N	-- : --
2	Lysis	1	710	00 : 05	fast	00 : 00	N	-- : --

Lysis temperature:65°C, lysis stop heating step: 2th

- 5) At the end of step 5, take the 96-deep-well pre-packed plate out. Add 300μL isopropanol into column 1 and 7, and put the 96-deep-well plate back into NPA-32p nucleic acid purification system. The automatic extraction experiment was carried out according to the following procedure:

Step	Name	Well	Volume (μL)	Waiting Time (min : ss)	Speed	Mixing Time (min : ss)	Adsorption	Magnet Time (min : ss)
1	Beads	6	200	-- : --	slow	00 : 30	N	00 : 30
2	Bind	1	1000	-- : --	fast	08 : 00	S	02 : 00
3	Wash 1	2	800	-- : --	fast	03 : 00	S	01 : 00
4	Wash 2	3	600	-- : --	fast	03 : 00	S	01 : 00
5	Wash 3	4	600	-- : --	fast	03 : 00	S	01 : 00
6	Elution	5	100	05 : 00	fast	10 : 00	C	01 : 40
7	Discard	6	200	-- : --	slow	00 : 30	N	-- : --

Elution temperature: 60°C, Elution start heating step: 6th

3. Oral swab processing

- 1) Use a cotton swab in intraoral wipes 15-20 times. Put the swab to a 2mL centrifuge tube, Cut out stem part from the swab with scissors.
- 2) Add 400μL—500μL of TES Buffer and 10 μL PK Solution to the 1.5mL DNase-free tube from the above.
- 3) Mix by pulse-vortexing for 15 seconds. Put it in a 56°C thermostat oscillator, Vortex 900rpm for 30mins. Centrifuge for a few seconds.
- 4) Transfer all liquid to the column 1and 7 of 96-deep-well plate.
- 5) Place 96-deep-well plate to the instrument, then plug in 8-strip tips and run the program.

Step	Name	Well	Volume (μL)	Waiting Time (min : ss)	Speed	Mixing Time (min : ss)	Adsorption	Magnet Time (min : ss)
1	Lysis	1	710	-- : --	fast	10 : 00	N	-- : --
2	Lysis	1	710	00 : 05	fast	00 : 00	N	-- : --

Lysis temperature:65°C, lysis stop heating step: 2th

- 6) At the end of step 5, take the 96-deep-well pre-packed plate out. Add 300μL isopropanol into column 1 and 7, and put the 96-deep-well plate back into NPA-32P nucleic acid purification system. The automatic extraction experiment was carried out according to the following procedure:

Step	Name	Well	Volume (μL)	Waiting Time (min : ss)	Speed	Mixing Time (min : ss)	Adsorption	Magnet Time (min : ss)
1	Beads	6	200	-- : --	slow	00 : 30	N	00 : 30
2	Bind	1	1000	-- : --	fast	08 : 00	S	02 : 00
3	Wash 1	2	800	-- : --	fast	03 : 00	S	01 : 00
4	Wash 2	3	600	-- : --	fast	03 : 00	S	01 : 00
5	Wash 3	4	600	-- : --	fast	03 : 00	S	01 : 00
6	Elution	5	100	05 : 00	fast	10 : 00	C	01 : 40
7	Discard	6	200	-- : --	slow	00 : 30	N	-- : --

Elution temperature: 60°C, Elution start heating step: 6th

4. Blood spot processing

- 1) Put 8 pieces of blood spots sample to a 1.5mL DNase-free tube. The area of each blood spots is 3×3mm.
- 2) Add 300μL TES Buffer and 10μL PK into the 1.5mL DNase-free tube.
- 3) Mix by pulse-vortexing intensively for 10 seconds. Put it in a 56°C thermostat oscillator, Vortex 900rpm for 45mins.
- 4) Transfer all liquid to the column 1 and 7 of 96-deep-well plate. (**Discard 150μL Lysis buffer from column**