



FastGene® 50 bp DNA Marker

Cat. No.	Product	Size	Concentration
MWD50	FastGene® 50 bp DNA Marker	100 reactions	74.05 µg / 500 µl

Storage conditions

Stable for up to 12 months at 4°C.

Stable for up to 24 months at -20°C.

Description

FastGene® 50 bp DNA Marker is a 17 fragments molecular weight standard for agarose gel electrophoresis. Using the unique combination of a number of proprietary plasmids digested with appropriate restriction enzymes and PCR products yields fragments ranging from 50-1,500 base pairs. The 4 bands at 500, 1000, 1200, and 1500 base pairs have specifically increased in intensity to serve as the reference bands. The approximate mass of DNA in each band is provided (0.7405 µg per load) for approximating the mass of DNA in comparably intense samples of similar size.

Kit content(s)

FastGene® 50 bp DNA Marker	500 µl x 1 vial
----------------------------	-----------------

Source

PCR products and double-stranded DNA digested with appropriate restriction enzymes are phenol extracted and equilibrated to 10 mM Tris-HCl (pH 8.0) and 10 mM EDTA.

Required materials but not provided

- Horizontal Electrophoresis system
- Power supplies
- Vortex or equivalent
- Microcentrifuge

Reaction Setup

5 µl/well, containing bromophenol blue & xylene cyanol FF as tracking dyes.


Range

50 - 1,500 bp



Number of bands

17



2% TAE agarose gel

NIPPON Genetics EUROPE GmbH

Mariaweilerstraße 28 a, 52349 Düren, Germany

Telephone: +49 2421 55 496 0

Fax: +49 2421 55 496 11