



Instruction Manual ibidi Heating System, Multi-Well Plates, K-Frame

Version 2.0



10929 ibidi Heating System, Multi-Well Plates, K-Frame

Contact

ibidi GmbH

Lochhamer Schlag 11
82166 Gräfelfing
Germany

Phone: +49 89 / 520 46 17 - 0

Fax: +49 89 / 520 46 17 - 59

E-mail: info@ibidi.de

Internet: ibidi.com

Contents

1	Preamble	5
1.1	Introduction	5
1.2	Safety Symbols	5
1.3	Nomenclature	6
1.4	Specifications	6
1.5	Disclaimer	8
1.6	Safety Considerations	8
1.7	Regulatory Statement	9
1.8	Limited Warranty	10
1.9	Transporting the ibidi Heating System, Multi-Well Plates, K-Frame	10
1.10	Repairing the ibidi Heating System, Multi-Well Plates, K-Frame	11
1.11	Waste Disposal – WEEE/RoHS Compliance Statement	11
1.11.1	EU Directive WEEE	11
1.11.2	EU Directive RoHS	12
2	Intended Use	13
3	Principle	13
4	Equipment	14
4.1	Components of the ibidi Heating System, Multi-Well Plates, K-Frame	14
4.2	Combination Options	15
4.3	Temperature Controller	16
4.4	Heated Lid	17
4.5	Heated Plate	18

5	Operation	19
5.1	Installation and Connection of the Parts	19
5.2	Start the ibidi Heating System, Multi-Well Plates, K-Frame	20
5.3	Setting Parameters in the Front Display	20
5.3.1	Run/Stop Mode	21
5.3.2	Incubation Parameters	21
5.3.3	Setup Menu	21
5.4	Temperature Adjustment in the Sample	23
5.5	Sample Preparation	24
6	IncuControl Software	25
7	Maintenance	26
7.1	Disinfection and Cleaning	26
7.2	Influence of Ambient Temperature and Ventilation	26
8	Troubleshooting	27
8.1	Focus not Stable	27
8.2	Evaporation Is too High	27
8.3	Condensation Inside the Stage Top Incubator	27

1 Preamble

1.1 Introduction

This manual is your guide to using the ibidi Heating System, Multi-Well Plates, K-Frame for cell culture experiments on an optical microscope. It instructs first-time users how to use the instrument, and serves as a reference for experienced users.

Before using the ibidi Heating System, Multi-Well Plates, K-Frame, please read this instruction manual carefully and make sure that the contents are fully understood. This manual should be easily accessible to the operator at all times during instrument operation. If this manual gets lost, order a replacement from www.ibidi.com.

To ensure safe operation, the ibidi Heating System, Multi-Well Plates, K-Frame must only be operated with the supplied components and according to the instruction manual.

1.2 Safety Symbols

Note that the signal words **WARNING**, **CAUTION** and **NOTE** have specific meanings in this manual. Do not proceed beyond a signal word until you have performed the indicated actions.

WARNING! A potentially hazardous situation which, if not avoided, could result in serious injury or even death. Warning messages in the text are displayed in a gray shaded box.

CAUTION A potentially hazardous situation which, if not avoided, could result in minor or moderate injury. It is also used to alert against damaging the equipment or the instrument.

NOTE Additional information to help achieve optimal instrument and assay performance.

Symbols on the product identification label and back panel of the device:



CE Marking: This symbol indicates the product’s compliance with EU legislation.



This label is positioned on the back of the device and prompts you to read the manual before using the device.




Product disposal: The symbol indicates that this product must be recycled/disposed of separately from other household waste. See page 11 for details.


1.3 Nomenclature



Temperature Controller



Heated Lid



Heated Plate with Heated Glass Bottom

1.4 Specifications

Table 1 – Specifications of the ibidi Heating System, Multi-Well Plates, K-Frame

Electrical Specifications Power Supply	
Protection class	I
International protection marking (IEC 60529)	IP 20
Overvoltage category	II
External power supply	AC 100 – 240 V, 50/60 Hz, 2 A
Input line voltage Temperature Controller	DC 24 V, 6.67 A, 160 W
Output voltage to channel 1 (Heated Lid), K-Frame	DC 10 V, max. 5 A
Output voltage to channel 2 (Heated Plate), K-Frame	DC 24 V, max. 5 A
Output voltage to channel 3 (Heated Glass Bottom), K-Frame	DC 10 V, max. 5 A
Operating Conditions	
Operating site	Indoor use only
Operating temperature	18-30°C/64-86°F (min 5°C/9°F less than set temperature)
Humidity	80% relative humidity (RH) up to 31°C 30% relative humidity (RH) up to 40°C
Operating altitude	max. 2000 m (atmospheric pressure 800-1060 hPa/11.6-15.4 psi)
Storage Conditions	-5-50°C/23-122°F, humidity <60% relative humidity (RH)
Outer Dimensions and Characteristics of the Components	
Temperature Controller	90 mm × 170 mm × 230 mm

Table 1 – (continued)

Heated Lid	Weight: 1720 g/3.8 lbs 25 mm × 160 mm × 108 mm Length of cable: 1.5 m Connector to Gas Incubation: Female Luer Lock
Heated Plate	9 mm × 160 mm × 130 mm K-Frame fitting: 160 mm × 110 mm Observation area: 107.5 mm × 72.5 mm Length of cable: 1.5 m
Heated Plate assembled with Heated Lid	Height: 30 mm Weight: 430 g/0.95 lbs
USB cable	Length: 1.8 m
Power supply cable	Length: 2.0 m (power supply to wall) Length: 1.2 m (power supply to device)

Temperature Control Range

Heated Lid (Channel 1)	Ambient temperature (min. 18°C) to +45°C
Heated Plate (Channel 2)	Ambient temperature (min. 18°C) to +45°C
Heated Glass Bottom (Channel 3)	Ambient temperature (min. 18°C) to +45°C
Temperature deviation	max. 2°C between wells of a Well-Plate

Recommended Temperatures before Adjustment

See Section 5.4

Optical Properties Heated Lid and Heated Glass Bottom

Glass	Selected HQ Floatglass
Thickness of the glass plate	1.1 mm/0.7 mm
Refractive index glass	$n_D = 1.520$ (588 nm)
Thickness of ITO coating	100 nm
Passivation layer	Thickness: 20-25 nm
Refractive index ITO with passiv. layer	1.95

Requirements for Objective Lenses

Minimum 1 mm working distance for coverslip bottom plates, e.g., ibidi μ -Plate 96 Well
 Minimum 2 mm working distance for thick bottom plates
 The sample focal plane in the K-Frame lies up to 3.5 mm above the normal sample focal plane. Therefore the objective must have at least 3.5 mm of travel distance along the focal axis to be able to image the sample.

Microscope Requirements

Working distance condenser	≥ 26 mm
Stage opening	160 mm × 110 mm

1.5 Disclaimer

- ibidi shall not be held liable, either directly or indirectly, for any damage incurred as a result of product use.
- The contents of this manual are subject to change without notice for product improvement.
- This manual is considered complete and accurate at publication.
- This manual does not guarantee the validity of any patent rights or other rights.
- If an ibidi software program doesn't function properly, this may be caused by a conflict from another program operating on the computer. In this case, take corrective action by uninstalling the conflicting product(s).
- ibidi is a registered trademark of ibidi GmbH in Germany and other countries.

1.6 Safety Considerations

WARNING!

- Only operate the ibidi Heating System, Multi-Well Plates, K-Frame with the supplied components.
- Only use the cables and plugs delivered with the system. The power plug of the control unit must be inserted in an outlet with a ground (earth) contact.
- Do not replace detachable power cables by power cables with inadequate specifications. By violating these instructions you risk electric shock and fire.
- Only use extension cables that have a protective ground wire.
- Do not operate the ibidi Heating System, Multi-Well Plates, K-Frame under conditions that pose a risk of explosion, implosion, or the release of gases. Only operate the ibidi Heating System, Multi-Well Plates, K-Frame with aqueous solutions.
- Do not operate a damaged ibidi Heating System, Multi-Well Plates, K-Frame. If the housing seems damaged or something is rattling inside the controller, contact the [ibidi service hotline](#) for repair.
- Some accessible parts of the Heated Plate and Heated Lid can reach temperatures up to 55°C. Avoid touching the temperature-controlled parts of the system when you have set the Temperature Controller to high temperatures.

CAUTION

- Ensure that the external power supply is easily accessible. The ibidi Heating System, Multi-Well Plates, K-Frame must be installed in a manner such that none of its components hinders access to the external power supply.

- Immediately replace damaged cords, plugs, or cables to avoid risk of personal injury or damage to the instrument.
- Only ibidi technical staff and technical staff instructed by ibidi are permitted to open and service the ibidi Heating System, Multi-Well Plates, K-Frame.
- The external power supply should not be brought into contact with moisture. If the housing is damaged, the external power supply should not be used.
- Avoid strong magnetic fields and sources of high frequency. The ibidi Heating System, Multi-Well Plates, K-Frame might not function properly when located near a strong magnetic field or high frequency source.
- Avoid vibrations from vacuum pumps, centrifuges, electric motors, processing equipment, and machine tools.
- Avoid dust and corrosive gas. Do not install the ibidi Heating System, Multi-Well Plates, K-Frame where it could be exposed to high levels of dust or to outside air or ventilation outlets.
- Install the ibidi Heating System, Multi-Well Plates, K-Frame in a horizontal and stable position, such as a table, bench, or desk upon which the instrument is installed.
- Install the ibidi Heating System, Multi-Well Plates, K-Frame in a location that enables easy access for maintenance.
- Do not place heavy objects on the instrument.
- The heated glass plates of the incubation chamber can break on mechanical impact. If so, the glass shards can lead to injuries if handled.
- Be aware that when switched on, a 10 V DC voltage is applied to the underside of the glass on the Heated Lid. Do not touch the underside or put it in contact with anything conductive. This could cause a short circuit that may damage the Temperature Controller and/or the Heating Devices.

1.7 Regulatory Statement

The ibidi Heating System, Multi-Well Plates, K-Frame has been designed, produced and tested in compliance with the European standard DIN EN 61010-1 (IEC 61010-1, "Safety requirements for electrical equipment for measurement, control and laboratory use"). Furthermore it meets the IEC 61326-1 ("Electrical equipment for measurement, control and laboratory use - EMC requirements") and CISPR 11 ("International Standard for electromagnetic emissions (disturbances) from Industrial, Scientific and Medical (ISM) Equipment") standards .

The device carries the CE mark.

The ibidi Heating System, Multi-Well Plates, K-Frame meets the Low Voltage Directive 2014/35/EU and the EMC Directive 2014/30/EC.

1.8 Limited Warranty

Products manufactured by ibidi, unless otherwise specified, are warrantied for a period of one year from the date of shipment to be free of defects in materials and workmanship. If any defects in the product are found during this warranty period, ibidi will repair or replace the defective part(s) or product free of charge.

This warranty does not apply to defects resulting from the following:

1. Improper or inadequate installation.
2. Improper or inadequate operation, maintenance, adjustment, or calibration.
3. Unauthorized modification or misuse.
4. Use of unauthorized tubing or fluidic connectors.
5. Use of consumables, disposables, and parts not supplied by an authorized ibidi distributor.
6. Corrosion due to the use of improper solvents, samples, or due to surrounding gases.
7. Accidents beyond ibidi's control, including natural disasters.

This warranty does not cover consumables, such as cell culture chambers and dishes, tubes, fluidic connectors, reagents etc.

The warranty for all parts supplied and repairs provided under this warranty expires on the warranty expiration date of the original product.

1.9 Transporting the ibidi Heating System, Multi-Well Plates, K-Frame

The weight of the Temperature Controller is approx. 1.7 kg/3.8 lbs. Moving the Temperature Controller during operation can pose a risk of personal injury or damage to the instrument.

For transport, switch off the Temperature Controller and then disconnect the heated components from the controller. Carry the devices carefully and avoid mechanical shocks.

WARNING!

Hot surface (max. 55°C)! Do not touch Heated Lid and Heated Plate when hot. Always disconnect the instrument from the power supply before transport and leave the instrument to cool down for approx. 5 minutes.

1.10 Repairing the ibidi Heating System, Multi-Well Plates, K-Frame

For inquiries concerning repair service, contact the ibidi service personnel and provide the model name and serial number of your system.

ibidi GmbH
Service Hotline: service@ibidi.com

CAUTION

Do not try to repair the ibidi Heating System, Multi-Well Plates, K-Frame by yourself. Disassembly of the ibidi Heating System, Multi-Well Plates, K-Frame is not allowed. Disassembly poses a risk of personal injury or damage to the devices. Contact ibidi service personnel if there is a need to disassemble a device.

1.11 Waste Disposal – WEEE/RoHS Compliance Statement

The European Union (EU) has enacted two directives, the first on product recycling (Waste Electrical and Electronic Equipment, WEEE) and the second on limiting the use of certain substances (Restriction on the use of Hazardous Substances, RoHS).

1.11.1 EU Directive WEEE

The ibidi Heating System, Multi-Well Plates, K-Frame must be disposed of in compliance with the WEEE Directive 2012/19/EC.



This symbol on the product is in accordance with the European Union’s Waste Electrical and Electronic Equipment (WEEE) Directive. The symbol indicates that this product must be recycled/disposed of separately from other household waste. It is the end user’s responsibility to dispose of this product by taking it to a designated WEEE collection facility for the proper collection and recycling of the waste equipment. The separate collection and recycling of waste equipment will help to conserve natural resources and protect human health and the environment. For more information about recycling, please contact your local environmental office, an electrical/electronic waste disposal company or distributor where you purchased the product.

1.11.2 EU Directive RoHS

Two Categories of products covered by the WEEE Directive are currently exempt from the RoHS Directive – Category 8, medical devices (with the exception of implanted or infected products) and Category 9, monitoring and control instruments.

All of our products fall into either Category 8 or 9, and are currently exempt from the RoHS Directive. Nevertheless, the ibidi Heating System, Multi-Well Plates, K-Frame meets the requirements set forth in the RoHS Directive 2011/65/EC.

2 Intended Use

The ibidi Heating System, Multi-Well Plates, K-Frame is a stage top incubator for live cell imaging, that fits onto microscope stages with a mounting frame K opening (110 × 160 mm). The geometry of the chamber with a heated glass plate underneath and also above the sample, provides a platform for keeping multi-well plates at a constant temperature during the experiment on the microscope.

For full control of the incubation environment, the ibidi Heating System, Multi-Well Plates, K-Frame can be combined with the ibidi Gas Incubation System to regulate the CO₂ and O₂ concentrations as well as the humidity.

3 Principle

The Heated Lid and Heated Plate with the Heated Glass Bottom are designed to keep cells on-stage in a multiwell plate at 37°C. The Heated Plate is mounted on the microscope stage. The multiwell plate is placed directly on the Heated Glass Bottom and the Heated Lid is put on top.

The gas mixture is streamed into the chamber through the gas inlet and then slowly escapes through the gaps between the Heated Lid, Heated Plate, and Heated Glass Bottom (small gray arrows in Figure 1).




Figure 1 – Schematic cross-sectional view of the ibidi Heating System, Multi-Well Plates, K-Frame.

The Heated Lid prevents condensation effects inside the entire incubation system (Figure 2).




Figure 2 – The Heated Lid prevents condensation effects on the lid of the cell culture vessel.








4 Equipment

The ibidi Heating System, Multi-Well Plates, K-Frame consists of the Heated Lid, Heated Plate with Heated Glass Bottom, and a Temperature Controller.

This section gives an overview of the components and their functionality and characteristics.

4.1 Components of the ibidi Heating System, Multi-Well Plates, K-Frame

The components of the ibidi Heating System, Multi-Well Plates, K-Frame (#10929) are listed below.

Description	Drawing
Temperature Controller	
Heated Lid with electrical cable and D-sub connector to connect to the Temperature Controller	
Heated Plate with Heated Glass Bottom with electrical cable and D-sub connector to connect to the Temperature Controller	
Temperature Adjustment Set: 1 temperature sensor 1 perforated ibidi μ -Plate 24 Well	
USB cable to connect the Temperature Controller with a computer Country specific power cord to connect the external power supply to the wall socket	
External power supply for the Temperature Controller	
USB flash drive with IncuControl software and instruction manual	

4.2 Combination Options

The parts of the ibidi Heating System, Multi-Well Plates, K-Frame are combined as shown in Figure 3. The Heated Lid and Heated Plate fit on an inverted microscope stage equipped with a mounting frame K (160 × 110 mm).

The Heated Plate can hold all well-plates matching the ANSI SLAS 1-2004 standard (R2012) with a maximum height of 23 mm.

The ibidi Heating System, Multi-Well Plates, K-Frame can be combined with the Gas Incubation System, that provides CO₂ and O₂ (optional) control and a defined humidified atmosphere. Detailed information on the Gas mixer unit is given in the Gas Incubation System instructions.




Figure 3 – Overview of the parts of the ibidi Heating System, Multi-Well Plates, K-Frame with combination options

4.3 Temperature Controller

The Temperature Controller is designed to control the different heated components (four channels available).

The front display shows the set values (S) and the current values (I). The settings can be adjusted via the control buttons (Section 5.3).

If you wish to control the Temperature Controller via PC, use the IncuControl software (Section 6).




Figure 4 – Front view of the Temperature Controller.

The LEDs on the front indicate the status of the channels, connection to the power supply, and USB connection (Figure 4).

Control LEDs Channel 1-4	
LED off	Channel inactive
LED on	Channel active
LED fast blinking	Channel error

All plugs for the electrical connections are integrated into the rear of the Temperature Controller (Figure 5). The setup of the connections is shown in Section 9.




Figure 5 – Rear view of the Temperature Controller.

4.4 Heated Lid

The Heated Lid provides excellent optical quality, and also allows for the use of all standard microscopy techniques, including differential interference contrast (DIC). Due to the height of the lid, we recommend using condensers with a working distance of ≥ 26 mm.

The upper glass part of the Heated Lid is electrically heat controlled. The Heated Lid fits exactly and securely onto the Heated Plate.

For gas incubation, the Heated Lid is equipped with inlets for the gas flow and the humidity sensor. Detailed information on the Gas Incubation System is given in the separate instructions.




Figure 6 – Parts of the Heated Lid

4.5 Heated Plate

The Heated Plate provides the base for the ibidi Heating System, Multi-Well Plates, K-Frame. The lower part of the Heated Plate fits into a mechanical or scanning stage with opening for mounting frames K (160 mm × 110 mm).




Figure 7 – Heated Plate, top view.

The Heated Plate can be leveled to compensate uneven samples. Use the 0.9 mm INBUS/Hex Key to adjust the M2 set screws at each of the plate's corners.




Figure 8 – Leveling of the Heated Plate.

5 Operation

Before starting an experiment, check that the ibidi Heating System, Multi-Well Plates, K-Frame fits on your microscope stage, and that your multi-well plates are compatible to the Heated Plate. Connect all parts (Section 5.1) and perform a temperature adjustment as explained in Section 5.4.

Important!

The Temperature Controller only controls the temperature at the internal sensors placed on the Heated Lid and in the Heated Plate. The resulting temperature in the sample must be adjusted for your specific setup. Follow the instructions in Section 5.4.

5.1 Installation and Connection of the Parts

The components of the ibidi Heating System, Multi-Well Plates, K-Frame are connected as shown in Figure 9.




Figure 9 – Installation of the components of the ibidi Heating System, Multi-Well Plates, K-Frame.

1. Place the Temperature Controller next to the microscope and connect the power supply.
2. Confirm that the power switch of the Temperature Controller is off.
3. Insert the Heated Plate into the opening of the microscope stage.
4. Place the Heated Lid onto the Heated Plate.

5. Connect the electrical cables of the Heated Lid (Channel 1), Heated Plate (Channel 2), and Heated Glass Bottom (Channel 3) to the plugs on the back of the Temperature Controller.

To adjust the sample temperature, use the temperature sensor of the Temperature Adjustment Set and plug it into the corresponding connector on the back of the Temperature Controller. The other end is placed in the sample (Section 5.4).

Optional: To setup communication with the IncuControl Software, the USB cable must be connected between the Temperature Controller and the computer.

5.2 Start the ibidi Heating System, Multi-Well Plates, K-Frame

The Temperature Controller is switched on by the dip-switch on the back. Make sure the display shows "RUN" (see Section 4.3). The system immediately starts heating up the heated components.

Warning!

Be aware that when the system is switched on, 10V DC voltage is applied to the underside of the glass plates. Do not touch the underside or contact it with anything conductive! This could cause a short circuit that may destroy the controller and/or the lid.

If you are not sure which temperature is set (e.g. when operating the system the first time), it is recommended to disconnect all heated devices (Heated Plate, Heated Lid, and other optional heated devices) and then switch the system on. By doing this, you ensure that you do not start the heating process with the wrong temperature settings. It is now possible to set the temperatures for the individual channels (Section 5.3). After this, you can re-connect the heated devices.

5.3 Setting Parameters in the Front Display

All control parameters can be manually set on the controller using the buttons and the display on the front panel.




Figure 10 – Temperature Controller front display and set buttons.

The display shows the measured ('I' = instantaneous) and set ('S' = set) temperatures of all the channels.

The cursor position is indicated with square brackets ("[]"). You can move the cursor using the "left" and "right" buttons. If you want to select a parameter or a function, press the round button and the square brackets will change to angle brackets ("< >"). Now you are able to change the value with the "up" and "down" buttons. To confirm the changed value, you must press the round button once more.

Standardly, the channels are pre-defined with the following settings to match all possible combinations of heated components on the controller:

Table 4 – Channel assignment of the Temperature Controller

Channel No.	Short name	Components to connect	Product number
Channel 1	"Lid"	Heated Lid, Universal Fit or Heated Lid, Multi-Well Plates, K-Frame	#10918, #10927 #10929
Channel 2	"Plate"	Heated Plate, Universal Fit, for 1 Chamber or Heated Plate, Multi-Well Plates, K-Frame	#10918 #10929
Channel 3	"Glass"	Heated Glass Bottom, Multi-Well Plates, K-Frame	#10929
Channel 4	"Ch4"	Heated Plate, Universal Fit, for 4 µ-Slides	#10927, #10928

5.3.1 Run/Stop Mode

Set the whole system to run or stop mode by manipulating the setting in the display's bottom line to "Run" (= system is running) or "Stop" (= system is not running).




Figure 11 – Main display indicating "Run"

5.3.2 Incubation Parameters

Set the parameters in the front display and wait for the temperature to equilibrate. It is recommended to start the system at least 30 minutes before inserting cells. Recommended values are given in table 5.

Before starting an experiment, an initial temperature adjustment must be performed (see Section 5.4).

5.3.3 Setup Menu

The Setup menu offers the possibility to change the settings of the individual channels and the display. For standard operation no changes need to be made. Enter the setup menu by navigating to the "Setup" entry and pressing the round button. The mode dialog will open.

Mode Each channel of the Temperature Controller can be set to an "On" or "Off" state. For standard operation all channels must be turned on.




Figure 12 – Mode dialog

1. Press the round button once more to move the cursor to the channel number.
2. Select the respective channel and confirm with the round button.
3. Select "On" or "Off" and confirm with the round button.
4. To move on navigate to "Return" and press the round button.

Alarms In this dialog it is possible to set the alarm limits for divergent control parameters. Move from the mode menu to the alarm menu by pressing the "right" button.

The maximum and minimum limits of the alarm can be set for each individual channel. If the alarm is activated (which happens when the current value goes under/over the low/high limits), the current value (I) blinks, showing alternately the value and "low" or "high", respectively. To stop the display blinking, navigate to the blinking channel and press the round button once.




Figure 13 – Alarm settings dialog

Preferences Set the brightness and contrast of the display in the preferences menu.




Figure 14 – Preferences dialog

Info Info about serial number and firmware version is shown in the info dialog.




Figure 15 – Info display

5.4 Temperature Adjustment in the Sample

The Temperature Controller controls the temperature of the actively heated components. The temperature in the sample is affected primarily by the temperature of the Heated Plate and the Heated Glass Bottom, and, due to heat loss to the surrounding environment, they must usually be set to a temperature slightly above the desired temperature in the sample. The Heated Lid acts primarily to prevent greater heat loss to the environment and to prevent condensation in the chamber during microscopy experiments; it should generally have a temperature several degrees above that of the Heated Plate. To ensure that the sample reaches the correct temperature, ibidi recommends a temperature adjustment step before starting the first experiment. The temperature adjustment should also be done separately for each chamber type used.

When starting the temperature adjustment in the sample in your specific setup, use the recommended temperature settings:

Table 5 – Recommended operation parameters

Parameter	Recommended Values
Heated Lid	42°C
Heated Plate	38°C
Heated Glass Bottom	37°C

To measure and calibrate the sample temperature, use the provided Temperature Adjustment Set. The temperature of the heated components should be adjusted in small steps. Follow this procedure to adjust the sample temperature:

1. Connect the temperature sensor (thermocouple type K) to the plug labeled "Temp Sensor" on the Temperature Controller's back.
2. Fill the provided multi-well plate with water. The volume of water and type of plate should correspond to those of your experiment.
3. Put the loose end of the thermocouple through the hole in the lid of the multi-well plate, making sure that it is submerged in water.
4. Place the multi-well plate onto the Heated Plate.
5. Set all external parameters, such as room temperature, air conditioning, airflow, illumination, microscope settings, etc, to the values that you will use in your experiment. For a 37°C sample temperature, use the recommended temperatures in table 5.3.2.
6. After 30 minutes, check the temperature of the external sensor on the display or in the IncoControl software. If the sample temperature is still too low, raise the set values for the Heated Lid, Heated Plate, and Heated Glass Bottom all uniformly by 0.5-1°C and wait until the sample temperature is stable again (at least 10-15 min).
7. When the sample temperature has reached the desired value, write down set values for Heated Plate, Heated Glass Bottom, and Heated Lid.

Note!

The temperature of the Heated Lid must be set several degrees warmer than the temperature of the Heated Plate!

When changing the incubation temperature, the Heated Glass Bottom must always be set to 1°C lower than the temperature of the Heated Plate, e.g. if the Heated Plate is set to 39°C, the Heated Glass Bottom must be set to 38°C. If the Heated Glass Bottom isn't set 1°C lower, the temperature distribution in the sample may not achieved the specified values.

The temperature adjustment must be repeated from time to time (at least once a year), especially if one of the following conditions has been changed:

- Room temperature
- Air conditioning
- Chamber type or objective lens
- Humidification and gas flow
- Use of an XL-Incubator
- Use of an Objective Heater

5.5 Sample Preparation

- Prepare the cells according to your protocol and place the multi-well plate in the Heated Plate mounted on the microscope stage.
- Close the Heated Lid.
- Let the temperature of the system equilibrate for a minimum of 30 min before you start your experiments.

Note!

Working with high aperture objectives might limit the observation area due to their large diameter. In this case ibidi recommends leaving out the outermost rows of wells to avoid the risk of objective collision with the Heated Plate.

6 IncuControl Software

The Temperature Controller has a USB interface for computer control and data logging. For this purpose, ibidi provides the IncuControl software that comes with the controller or can be downloaded from the [ibidi website](#).

For more details, please refer to the [IncuControl instructions](#).




Figure 16 – Temperature control window in IncuControl.

Prohibitions on the use of ibidi software:

- Copying software for other than backup purposes
- Transferring or licensing of the right to use software to a third party
- Disclosure of confidential information regarding software
- Modification of software

7 Maintenance

7.1 Disinfection and Cleaning

The heated components and the Temperature Controller do not have to be sterile, because the incubation chamber has no direct contact with the cells and the cell culture medium. If disinfection is necessary for some reason, we recommend using isopropanol (70%) or common lab disinfection solutions based on quaternary ammonium compounds (e.g. Barrycidal 36 or Pharmacidal).

All parts of the ibidi Heating System, Multi-Well Plates, K-Frame can be cleaned from the outside. We recommend using ultrapure water for cleaning. Fingerprints on the Heated Lid can be removed using isopropanol (70%) or lens/eyeglass cleaning wipes.

CAUTION

When cleaning the heated glass plate(s), be careful when wiping the inner surface of the lid not to damage the electro-conductive coating. Also take care not to damage the glass plate. The use of ethanol or other types of organic solvents may remove the instrument's paint.

To clean the heated components and/or the Temperature Controller switch off the Temperature Controller and disconnect all electrical cables. Leave the instrument to cool down for approx. 5 minutes.

WARNING!

Hot surface (max. 55°C)! Do not touch Heated Lid, Heated Plate, and Heated Glass Bottom when hot.

7.2 Influence of Ambient Temperature and Ventilation

The ambient temperature affects the conditions inside the incubation chamber. Devices, such as computers and camera controllers, can significantly heat up small rooms. In this case, we recommend equilibrating the room temperature to the typical experimental conditions at least 2-3 hours before starting the experiment.

Ventilation can enhance the effect of temperature and humidity changes in the vicinity of the incubation chamber. In a case where the airflow (e.g., air conditioning) cannot be stopped, we recommend protecting the microscope as much as possible.

The use of an XL-Incubator and/or an objective heater minimizes those effects and helps significantly to stabilize surrounding conditions.

8 Troubleshooting

8.1 Focus not Stable

Focus drift is detrimental for most microscopy experiments, especially long duration time-lapse experiments. Focus stability is mainly influenced by mechanical changes and temperature variations. Follow these recommendations to keep your cells in focus:

- Switch on all components (e.g., heating, gas incubation, computer, or other equipment) at least 60 minutes before starting the experiment.
- After you put the sample onto the microscope, wait 30 minutes before starting a time-lapse experiment to achieve temperature and immersion oil equilibration¹.
- Keep the room temperature as stable as possible. Air conditioning should either be working continuously or switched off.
- Do not change the temperature during the experiments. Avoid door/window openings, as this could rapidly change the temperature.
- Eliminate all sources of mechanical vibrations. Use a damped table for your microscope.

8.2 Evaporation Is too High

Depending on the incubating conditions, small volumes might evaporate quickly, especially during long-term experiments. If you have an actively controlled humidifying device (e.g. ibidi Gas Incubation System), increase the set value for relative humidity. Additionally, we suggest using silicone oil (e.g. Anti-Evaporation Oil, ibidi, 50051) to decrease evaporation.

Covering the medium with sterile silicone oil prevents all evaporation effects and is compatible with cell culture. Please do not use mineral oil, as this can be harmful to your cultureware.

Equilibrate oil and medium inside the incubator overnight. This step helps to avoid the formation of air bubbles, and pre-warms the solutions to 37°C. Afterwards, fill your slide with cells and medium. Cover the medium's surface with an appropriate amount of silicone oil. Don't drip the oil directly onto the surface, but let it run down the edges of the cell culture vessel. Details about avoiding evaporation are given on the ibidi web site in [Application Note 12 "Avoiding Evaporation"](#).

8.3 Condensation Inside the Stage Top Incubator

Check the temperature of the chamber (Heated Lid and Heated Plate). Make sure the humidity sensor is not in contact with the Heated Lid. In case of condensation, decrease the humidity and air-dry the incubator if necessary.

Please contact ibidi at techsupport@ibidi.com for further troubleshooting help.

¹If the experiment needs to be started immediately, either after placing the sample on the microscope or after closing the lid, we recommend checking the focus for 20 minutes. In the first few minutes after starting the experiment, temperature equilibration might influence the focus/z-position of the cells.



Certified ISO 9001:2008, EN ISO 13485:2012

ibidi GmbH

Lochamer Schlag 11
82166 Gräfelfing (München)
Germany

Toll free within Germany:
Phone: 0800 / 00 11 11 28
Fax: 0800 / 00 11 11 29

International calls:
Phone: +49 89 / 520 46 17 - 0
Fax: +49 89 / 520 46 17 - 59

E-Mail: info@ibidi.com
www.ibidi.com