# sticky-Slide Tissue



Instruction Manual



The sticky-Slide Tissue supports immunofluorescence and immunohistochemistry staining protocols of tissue sections on microscope glass slides.

This document applies to the following product:

80518 sticky-Slide Tissue

#### Material

The material of sticky-Slides is identical to that of  $\mu$ -Slides. All sticky-Slides are delivered sterilized and individually packed. Please keep in mind that sterility is lost when non-sterile substrates are used. The sticky-Slides are not autoclavable, as they are only temperature-stable up to 60 °C/140 °F.

The sticky bottom itself is a  $85\,\mu$ m biocompatible double-faced adhesive tape. The tape is covered by a protection film, which must be removed before usage.

### **Shipping and Storage**

The sticky-Slides are sterilized and sealed in a gas-permeable packaging. The shelf life under proper storage conditions (in a dry place, no direct sunlight) is outlined in the following table.

Conditions		
Shipping conditions	Ambient	
Storage conditions	RT (15–25 <i>°</i> C)	
Shelf Life		
sticky-Slides	36 months	

Geometry	

Please note that the outer dimensions of the sticky-Slide Tissue are smaller compared to the standard ISO 8037/1 slide format.

Specifications		
Outer dimensions	24.0 × 55.0 mm <sup>2</sup>	
Total height	7.8 mm	
Channel height	235 µm	
Channel volume	160 µl	
Ceiling thickness	1.67 mm	
Adapters	Female Luer	
Bottom	None	



The sticky-Slide Tissue is made for tissue sections smaller than approx. 10 mm x 30 mm or 15 mm x 20 mm, respectively.

# Surface Compatibility

The sticky-Slide Tissue is compatible with flat, clean, dust-free, fat-free surfaces, such as microscope slides and coverslips. Best results are achieved with completely dry surfaces. Please test your specific surface with a free sample from ibidi.com.

### Handling and Assembly

Follow the protocol below to mount the sticky-Slide Tissue to a microscope slide bearing a tissue section after deparaffinization and rehydration.

- 1. Deparaffinize and rehydrate your tissue.
- Make sure the size of the tissue section is smaller than approx. 10 mm × 30 mm or 15 mm × 20 mm, respectively.
- 3. Clean and dry the surface around the tissue section. Do not touch the tissue section in order to keep it flat and even.



 Remove the protection film from the sticky-Slide. First, scrape off a corner of the protection film using a scalpel or pointy tweezers. After that, grab this corner and remove the protection film.



5. Mount microscope slide and sticky-Slide by pressing firmly with your fingers (use gloves) until both parts are completely sealed. Make sure the tissue section is aligned with the channel. For a better visualization of the tissue section, a dark background is beneficial.



6. Fill the channel using a syringe or a pipet. Incline the slide and inject the liquid in an upwards motion so that air can escape above the liquid level. Optionally, connect Luer adapters and tubing for automated liquid exchange.





**TIP** – When filling the channel with a pipet, point the pipet tip directly onto the channel opening and inject the liquid.

**TIP** – Inclining the slide helps filling the channel without entrapped air bubbles. Keep in mind that air is much lighter than water. Let the air escape above the liquid. **TIP** – In case air bubbles become a problem during filling, use the final re-hydration solution (typically approx. 40–70% ethanol) to fill the channel. Subsequently, wash with water or buffer to remove the ethanol.



For liquid exchange, a continuous exchange is recommended without removing any liquid from the channel itself: Add new liquid from one side removing it from the other. For a thorough wash, flush the channel with approx. 3× the channel volume, i.e. 480 µl.



- Follow your staining protocol.
- Optionally, remove the sticky-Slide Tissue with a scalpel and continue with mounting for microscopy.



**TIP** – Instead of using a microscope glass slide, a coverslip may be used to support high resolution microscopy.

# Chemical Compatibility

The following table provides basic information on the chemical and solvent compatibility of the sticky-Slide Tissue. For a full list of compatible solvents and more information on chemical compatibility, visit ibidi.com/chemicals.

<b>Chemical / Solvent</b>	Compatibility
Methanol	Yes
Ethanol	Yes
Formaldehyde	Yes
Acetone	No
Mineral oil	Yes
Silicone oil	Yes
Immersion oil	See Section "Immer- sion Oil"

### **Immersion Oil**

The compatibility with immersion oil depends on the used substrate.

### For research use only!

Further information can be found at ibidi.com. For questions and suggestions, please contact us by e-mail at info@ibidi.com or by telephone at +49 (0)89/520 4617 0. © ibidi GmbH, Lochhamer Schlag 11, 82166 Gräfelfing, Germany.