# KAPA HyperPlus Library Preparation Kit 取扱説明書

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## KAPA HyperPlus Library Preparation Kit

KR1145 – v1.15

This Technical Data Sheet provides product information and a detailed protocol for the KAPA HyperPlus Library Preparation Kit.

### Contents

Product Description	2
Product Applications	2
Product Specifications	2
Shipping and Storage	2
Handling	2
Quality Control	2
Important Parameters	3
Automated Library Construction	3
Safe Stopping Points	3
Input DNA	4
Fragmentation Parameters	5
Adapter Design and Concentration	5
Post-ligation Cleanup	6
Paramagnetic SPRI® Beads and Reaction Cleanups	6
Size Selection	6
Library Amplification	7
Evaluating the Success of Library Construction	8
Process Workflow	11
Library Construction Protocol	12
Appendix 1: Size Selection	15
Appendix 2: Optimization of Fragmentation Parameters	16
Restrictions and Liabilities	18
Note to Purchaser: Limited Product Warranty	18
Note to Purchaser: Limited License	18

Kit Codes and Components		
<b>KK8510</b> <b>KK8511*</b> 8 libraries	KAPA Frag Enzyme	100 µL
	KAPA Frag Buffer (10X)	50 µL
	KAPA Frag Conditioning Solution	580 µL
	End Repair & A-Tailing Buffer	70 µL
	End Repair & A-Tailing Enzyme	30 µL
	Ligation Buffer	300 µL
	DNA Ligase	100 µL
	KAPA HiFi HotStart ReadyMix (2X)* Library Amplification Primer Mix (10X)*	250 µL 50 µL
<b>KK8512</b> <b>KK8513*</b> 24 libraries	KAPA Frag Enzyme	270 µL
	KAPA Frag Buffer (10X)	140 µL
	KAPA Frag Conditioning Solution	580 µL
	End Repair & A-Tailing Buffer	210 µL
	End Repair & A-Tailing Enzyme	90 µL
	Ligation Buffer	900 µL
	DNA Ligase	300 µL
	KAPA HiFi HotStart ReadyMix (2X)* Library Amplification Primer Mix (10X)*	690 µL 138 µL
<b>KK8514</b> <b>KK8515*</b> 96 libraries	KAPA Frag Enzyme	1.27 mL
	KAPA Frag Buffer (10X)	640 µL
	KAPA Frag Conditioning Solution	580 µL
	End Repair & A-Tailing Buffer	930 µL
	End Repair & A-Tailing Enzyme	400 µL
	Ligation Buffer	3.8 mL
	DNA Ligase	1.26 mL
	KAPA HiFi HotStart ReadyMix (2X)* Library Amplification Primer Mix (10X)*	3.0 mL 600 µL

\*KK8511, KK8513 and KK8515 are available for PCR-free workflows, and do not contain library amplification reagents.

### Quick Notes

- This kit provides a versatile, streamlined DNA fragmentation and library construction protocol. Libraries for Illumina® sequencing may be prepared from a wide range of DNA samples and inputs (1 ng – 1 µg) in 1.5 – 3 hours.
- The novel one-tube chemistry improves library yield and quality, particularly for FFPE and low-input libraries.
- The protocol is easy to automate. Generous reagent excesses are supplied in 96-reaction kits to accommodate automated liquid handling.
- Kits contain all the reagents for DNA fragmentation, library construction and high-efficiency, low-bias library amplification, except for adapters and SPRI beads. Kits without an amplification module are available for PCR-free workflows.
- The **Process Workflow** (p. 11) provides an overview of the DNA fragmentation and library construction process. **Appendix 1** provides a protocol for dual-SPRI size selection. **Appendix 2** contains guidelines for the optimization of fragmentation parameters.

## Product Description

The KAPA HyperPlus Library Preparation Kit provides a versatile, streamlined DNA fragmentation and library construction protocol for the rapid preparation of libraries for Illumina® sequencing. The novel one-tube protocol improves the efficiency and consistency of library construction, and yields libraries of similar or better quality than those produced with the KAPA Hyper Prep Kit from Covaris-sheared DNA. It outperforms tagmentation-based workflows in terms of robustness, flexibility and sequence coverage and uniformity.

The workflow combines enzymatic steps and employs minimal SPRI® bead cleanups, thereby reducing sample handling and overall library preparation time to 1.5 – 3 hours. The kit contains all of the enzymes and reaction buffers required for:

1. **enzymatic fragmentation**, which produces dsDNA fragments;
2. **end repair and A-tailing**, which produces end-repaired, 5'-phosphorylated, 3'-dA-tailed dsDNA fragments;
3. **adapter ligation**, during which dsDNA adapters with 3'-dTTP overhangs are ligated to 3'-dA-tailed molecules;
4. **library amplification (optional)**, which employs high-fidelity, low-bias PCR to amplify library fragments carrying appropriate adapter sequences on both ends.

The kit provides a single concentrated buffer and a single enzyme mixture for enzymatic fragmentation, as well as for each of the two library construction steps—but does not include adapters or SPRI beads required for cleanups after adapter ligation and library amplification.

In order to maximize sequence coverage uniformity, it is critical to minimize library amplification bias. KAPA HiFi DNA Polymerase is designed for low-bias, high-fidelity PCR, and is the reagent of choice for NGS library amplification.<sup>1,2,3,4</sup> KAPA HyperPlus Kits include KAPA HiFi HotStart ReadyMix (2X), a ready-to-use PCR mix comprising all the components for library amplification—except primers and template. Kits also include KAPA Library Amplification Primer Mix (10X), designed for the high-efficiency amplification of Illumina libraries flanked by the P5 and P7 flow cell sequences. Kits without the amplification module (KK8511, KK8513 and KK8515) are available for PCR-free workflows. They may also be combined with KAPA HiFi Real-Time Library Amplification Kits (KK2701 and KK2702), or with KAPA HiFi HotStart Uracil+ ReadyMix Kits (KK2801 and KK2802) for the amplification of bisulfite-converted libraries.

1. Oyola, S.O. et al. *BMC Genomics* 13, 1 (2012).
2. Quail, M.A. et al. *Nature Methods* 9, 10 – 11 (2012).
3. Quail, M.A. et al. *BMC Genomics* 13, 341 (2012).
4. Ross, M.G. et al. *Genome Biology* 14, R51 (2013).

## Product Applications

KAPA HyperPlus Library Preparation Kits are ideally suited for low- and high-throughput NGS library construction

workflows that require DNA fragmentation, end repair, A-tailing, adapter ligation and library amplification (optional). Kits are designed for library construction from a wide range of sample types and inputs (1 ng – 1 µg), and are compatible with complex, genomic DNA; low-complexity samples such as small viral genomes, plasmids, cDNA and long amplicons; and low-quality DNA such as FFPE samples.

The entire workflow is automation-friendly and may be incorporated into workflows for a wide range of NGS applications, including:

- Whole-genome shotgun sequencing
- Whole exome or targeted sequencing, using Roche Nimblegen™ SeqCap EZ, Agilent SureSelect, Illumina TruSeq™, or IDT xGen™ Lockdown™ Probes, or other hybridization capture systems
- RNA-seq (selected applications)
- Methyl-seq (in combination with KAPA HiFi HotStart Uracil+ ReadyMix for library amplification)

## Product Specifications

### Shipping and Storage

The enzymes provided in this kit are temperature sensitive, and appropriate care should be taken during shipping and storage. KAPA HyperPlus Library Preparation Kits are shipped on dry ice or ice packs, depending on the destination country. Upon receipt, immediately store enzymes and reaction buffers at -15°C to -25°C in a constant-temperature freezer. When stored under these conditions and handled correctly, the kit components will retain full activity until the expiry date indicated on the kit label.

### Handling

Always ensure that KAPA HyperPlus Library Preparation Kit components have been fully thawed and thoroughly mixed before use. The KAPA Hyper Prep End Repair & A-Tailing Buffer may contain a precipitate when thawed. Vortex this buffer thoroughly to ensure that it is fully homogenized before use. KAPA HiFi HotStart ReadyMix (2X) contains isostabilizers and may not freeze completely, even when stored at -15°C to -25°C. Nevertheless, always ensure that the ReadyMix is fully thawed and thoroughly mixed before use. Reaction master mixes prepared from the enzymes and buffers for fragmentation, end repair and A-tailing, as well as for ligation, are very viscous and require special attention during pipetting. Keep all enzyme components and master mixes on ice as long as possible during handling and preparation.

### Quality Control

All kit components are subjected to stringent functional quality control, are free of detectable contaminating exo- and endonuclease activities, and meet strict requirements with respect to DNA contamination. Please contact [support@kapabiosystems.com](mailto:support@kapabiosystems.com) for more information.

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## Process Workflow

See **Table 1** on p. 4 for recommended inputs of different types of DNA for different sequencing applications.


For samples containing EDTA: include the appropriate concentration of KAPA Frag Conditioning Solution (see **Table 2** on p. 4)

See **Table 3** on p. 5 for recommended adapter stock concentrations.

Size selection, employing SPRI® beads or an electrophoretic method, may be incorporated at this point, if appropriate.

**PCR-free workflows end here.**

Size selection, employing SPRI beads or an electrophoretic method, may be incorporated at this point, if appropriate.

## Recommended QC Metrics

- **Concentration and Q-ratio of input DNA** (*KAPA hgDNA Quantification and QC Kit*)  
for FFPE DNA only

- **Concentration of adapter-ligated libraries** (*KAPA Library Quantification Kit*)  
Adapter-ligated libraries diluted in the range of 1/500 to 1/10,000 usually fall within the dynamic range of the assay

- **Electrophoretic profile of amplified libraries**
- **Concentration of amplified libraries** (*KAPA Library Quantification Kit*)  
Amplified libraries diluted in the range of 1/10,000 to 1/200,000 usually fall within the dynamic range of the assay

## Library Construction Protocol

### Notes:

- Standard fragmentation parameters may not result in the optimal size distribution for libraries prepared from your specific DNA samples. For this reason, precious samples should not be used when evaluating this kit for the first time. Instead, parameters should be optimized with a non-precious, bulk DNA sample that is representative of the actual samples to be processed. Please refer to **Appendix 2: Optimization of Fragmentation Parameters** for guidelines.
- This protocol does not include size selection. Please refer to **Appendix 1: Size Selection** for a detailed dual-SPRI® size selection protocol that may be included after ligation or after amplification.

### 1. Enzymatic Fragmentation

1.1 Dilute the amount of dsDNA to be used for library construction as follows:

- In 10 mM Tris-HCl (pH 8.0 – 8.5), in a total of 35 µL if the DNA preparation does not contain EDTA.
- In the EDTA-containing buffer in which samples are currently suspended, in a total of 30 µL if the DNA preparation does contain EDTA.

1.2 If the DNA samples contain EDTA: prepare a sufficient volume of appropriately diluted KAPA Frag Conditioning Solution (5 µL per DNA sample plus excess). Refer to **Table 2** (p. 4) for guidelines on the dilution of the KAPA Frag Conditioning Solution.

To each reaction tube/well of a PCR plate with 30 µL of EDTA-containing DNA, add 5 µL of diluted KAPA Frag Conditioning Solution Mix by gentle vortexing or pipetting up and down.

1.3 Assemble each fragmentation reaction on ice by adding the rest of the components in the order shown below:

Component	Volume
Double-stranded DNA (with KAPA Frag Conditioning Solution, if needed)	35 µL
KAPA Frag Buffer (10X)*	5 µL
KAPA Frag Enzyme*	10 µL
<b>Total volume</b>	<b>50 µL</b>

\*The KAPA Frag Buffer and Enzyme may be pre-mixed and kept on ice prior to reaction setup, and dispensed as a single solution.

1.4 Vortex gently and spin down briefly. Return the plate/tube(s) to ice. Proceed immediately to the next step.

1.5 Incubate in a thermocycler, pre-cooled to 4°C and programmed as outlined below. A heated lid is not required for this step. If used, set the temperature of the heated lid to ≤50°C.

Step	Temp	Time
Pre-cool block	4°C	N/A
Fragmentation	37°C	See table below
HOLD	4°C	∞

Mode fragment length	Incubation time at 37°C*	Optimization range
600 bp	5 min	3 – 10 min
350 bp	10 min	5 – 20 min
200 bp	20 min	10 – 25 min
150 bp	30 min	20 – 40 min

\*These parameters are a good starting point for high-quality genomic DNA. Please refer to **Appendix 2: Optimization of Fragmentation Parameters** for guidelines on how to optimize fragmentation time and temperature, if needed.

1.6 Transfer reactions to ice, and proceed immediately to **Section 2**.

### 2. End Repair and A-tailing

2.1 In the same tube(s)/plate in which enzymatic fragmentation was performed, assemble each end repair and A-tailing reaction as follows:

Component	Volume
Fragmented, double-stranded DNA	50 µL
End Repair & A-Tailing Buffer*	7 µL
End Repair & A-Tailing Enzyme Mix*	3 µL
<b>Total volume</b>	<b>60 µL</b>

\*The buffer and enzyme mix should preferably be pre-mixed and added in a single pipetting step. Premixes are stable for ≤24 hours at room temperature, for ≤1 week at 4°C, and for ≤3 months at -20°C.

2.2 Vortex gently and spin down briefly. Return the reaction tube(s)/plate to ice. Proceed immediately to the next step.

2.3 Incubate in a thermocycler programmed as outlined below. A heated lid is required for this step. If possible, set the temperature of the heated lid to ~85°C (instead of the usual 105°C).

Step	Temp	Time
End repair and A-tailing	65°C*	30 min
HOLD	4°C**	∞

\*Both the fragmentation and end repair enzymes are inactivated at 65°C. When reactions are set up according to recommendations, additional fragmentation should be negligible. The brief period of end repair is sufficient for enzymatically fragmented DNA.

\*\*If proceeding to the adapter ligation reaction setup without any delay, the reaction may be cooled to 20°C instead of 4°C.

2.4 Proceed immediately to the next step.

### 3. Adapter Ligation

- 3.1 Dilute adapter stocks to the appropriate concentration, as outlined in **Table 3** on p. 5.
- 3.2 In the same plate/tube(s) in which end repair and A-tailing was performed, assemble each adapter ligation reaction as follows:

Component	Volume
End repair and A-tailing reaction product	60 µL
Adapter stock (concentration as required)	5 µL
PCR-grade water*	5 µL
Ligation Buffer*	30 µL
DNA Ligase*	10 µL
<b>Total volume</b>	<b>110 µL</b>

\*The water, buffer and ligase enzyme should preferably be premixed and added in a single pipetting step. Premixes are stable for ≤24 hours at room temperature, for ≤1 week at 4°C, and for ≤3 months at -20°C.

- 3.3 Mix thoroughly and centrifuge briefly.
- 3.4 Incubate at 20°C for 15 min.

**Note:** to achieve higher conversion rates and library yields, particularly for low-input samples, consider increasing the ligation time, to a maximum of 4 h at 20°C, or overnight at 4°C. Please note that longer ligation times may lead to increased levels of adapter-dimer. Adapter concentrations may have to be optimized if ligation times are extended significantly.

- 3.5 Proceed immediately to the next step.

### 4. Post-ligation Cleanup

- 4.1 In the same plate/tube(s), perform a 0.8X SPRI® cleanup by combining the following:

Component	Volume
Adapter ligation reaction product	110 µL
Agencourt® AMPure® XP reagent	88 µL
<b>Total volume</b>	<b>198 µL</b>

- 4.2 Mix thoroughly by vortexing and/or pipetting up and down multiple times.
- 4.3 Incubate the plate/tube(s) at room temperature for 5 – 15 min to bind DNA to the beads.
- 4.4 Place the plate/tube(s) on a magnet to capture the beads. Incubate until the liquid is clear.
- 4.5 Carefully remove and discard the supernatant.
- 4.6 Keeping the plate/tube(s) on the magnet, add 200 µL of 80% ethanol.
- 4.7 Incubate the plate/tube(s) on the magnet at room temperature for ≥30 sec.

- 4.8 Carefully remove and discard the ethanol.
- 4.9 Keeping the plate/tube(s) on the magnet, add 200 µL of 80% ethanol.
- 4.10 Incubate the plate/tube(s) on the magnet at room temperature for ≥30 sec.
- 4.11 Carefully remove and discard the ethanol. Try to remove all residual ethanol without disturbing the beads.
- 4.12 Dry the beads at room temperature for 3 – 5 min, or until all of the ethanol has evaporated. **Caution: over-drying the beads may result in dramatic yield loss.**
- 4.13 Remove the plate/tube(s) from the magnet.
- 4.14 Resuspend the beads:
  - in 25 µL of elution buffer (10 mM Tris-HCl, pH 8.0 – 8.5) to proceed with library amplification (Section 5), or
  - in 55 µL of elution buffer (10 mM Tris-HCl, pH 8.0 – 8.5) to proceed with dual-SPRI size selection (**Appendix 1**).
- 4.15 Incubate the plate/tube(s) at room temperature for 2 min to elute DNA off the beads.
- 4.16 Place the plate/tube(s) on a magnet to capture the beads. Incubate until the liquid is clear.
- 4.17 Transfer the clear supernatant to a new plate/tube(s) and proceed with library amplification or size selection, as appropriate. Store purified, adapter-ligated libraries at 4°C for 1 – 2 weeks, or at -20°C.

### 5. Library Amplification

**Note:** Please refer to **Important Parameters: Library Amplification** and Section 7 of the **KAPA NGS Library Preparation Technical Guide** for more information on optimizing library amplification.

- 5.1 Assemble each library amplification reaction as follows:

Component	Volume
2X KAPA HiFi HotStart ReadyMix	25 µL
10X KAPA Library Amplification Primer Mix*	5 µL
Adapter-ligated library	20 µL
<b>Total volume</b>	<b>50 µL</b>

\*Or other, suitable 10X library amplification primer mix. The recommended final concentration of each primer in the library amplification reaction is 0.5 – 4 µM. Also refer to **Important Parameters: Library Amplification**.

- 5.2 Mix thoroughly and centrifuge briefly.
- 5.3 Amplify using the cycling protocol outlined on the next page:

## 5. Library Amplification (continued)

Step	Temp	Duration	Cycles
Initial denaturation	98°C	45 sec	1
Denaturation	98°C	15 sec	Minimum number required for optimal amplification (Table 4 or 5)
Annealing*	60°C	30 sec	
Extension	72°C	30 sec	
Final extension	72°C	1 min	1
HOLD	4°C	∞	1

\*Optimization of the annealing temperature may be required for non-standard (i.e., other than Illumina® TruSeq™) adapter/primer combinations.

5.4 Store the plate/tube(s) at 4°C or -20°C for up to 72 hours, or proceed directly to **Step 6: Post-amplification Cleanup**.

## 6. Post-amplification Cleanup

6.1 In the library amplification plate/tube(s), perform a 1X SPRI® cleanup by combining the following:

Component	Volume
Library amplification reaction product	50 µL
Agencourt® AMPure® XP reagent	50 µL
<b>Total volume</b>	<b>100 µL</b>

6.2 Mix thoroughly by vortexing and/or pipetting up and down multiple times.

6.3 Incubate the plate/tube(s) at room temperature for 5 – 15 min to bind DNA to the beads.

6.4 Place the plate/tube(s) on a magnet to capture the beads. Incubate until the liquid is clear.

6.5 Carefully remove and discard the supernatant.

6.6 Keeping the plate/tube(s) on the magnet, add 200 µL of 80% ethanol.

6.7 Incubate the plate/tube(s) on the magnet at room temperature for ≥30 sec.

6.8 Carefully remove and discard the ethanol.

6.9 Keeping the plate/tube(s) on the magnet, add 200 µL of 80% ethanol.

6.10 Incubate the plate/tube(s) on the magnet at room temperature for ≥30 sec.

6.11 Carefully remove and discard the ethanol. Try to remove all residual ethanol without disturbing the beads.

6.12 Dry the beads at room temperature for 3 – 5 min, or until all of the ethanol has evaporated. **Caution: over-drying the beads may result in dramatic yield loss.**

6.13 Remove the plate/tube(s) from the magnet.

6.14 Thoroughly resuspend the beads in an appropriate volume of elution buffer (10 mM Tris-HCl, pH 8.0 – 8.5) or PCR-grade water. **Always use PCR-grade water if proceeding to target capture.**

**Note:** If proceeding with a second post-ligation cleanup, or dual-SPRI size selection (**Appendix 1**), resuspend the beads in 55 µL of elution buffer.

6.15 Incubate the plate/tube(s) at room temperature for 2 min to elute DNA off the beads.

6.16 Place the plate/tube(s) on a magnet to capture the beads. Incubate until the liquid is clear.

6.17 Transfer the clear supernatant to a new plate/tube(s) and proceed with size selection (refer to **Appendix 1**), library QC, target capture or sequencing, as appropriate. Store purified, amplified libraries at 4°C for 1 – 2 weeks, or at -20°C.

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## Restrictions and Liabilities

This technical data sheet is provided “as is” and Kapa Biosystems assumes no responsibility for any typographical, technical, or other inaccuracies. The document is subject to change, without notice, in future editions.

To the maximum extent permitted by applicable law, Kapa Biosystems disclaims all warranties, either express or implied, with regard to this technical data sheet and any information contained herein, including but not limited to the implied warranties of merchantability and fitness for a particular purpose. Kapa Biosystems shall not be liable for errors or for incidental or consequential damages in connection with the furnishing, use, or performance of this document or of any information contained herein.

This document might contain references to third party sources of information, hardware or software, products, or services and/or third party web sites (collectively the “Third-Party Information”). Kapa Biosystems does not control, and is not responsible for, any Third-Party Information. The inclusion of Third-Party Information in this document does not imply endorsement by Kapa Biosystems of the Third-Party Information or the third party in any way.

Kapa Biosystems is not responsible nor will be liable in any way for your use of any software or equipment that is not supplied by Kapa Biosystems in connection with your use of Kapa Biosystems products.

Kapa Biosystems does not in any way guarantee or represent that you will obtain satisfactory results from using Kapa Biosystems products as described herein. The only warranties provided to you are included in the Limited Warranty enclosed with this document. You assume all risk in connection with your use of Kapa Biosystems products.

## Note to Purchaser: Limited Product Warranty

Any product that does not meet the performance standards stated in the product specification sheet will be replaced at no charge. This warranty limits our liability to the replacement of the product. No other warranties of any kind, express or implied, including without limitation, implied warranties of merchantability or fitness for a particular purpose, are provided by Kapa Biosystems. Kapa Biosystems shall have no liability for any direct, indirect, consequential or incidental damages arising out of the use, the results of use or the inability to use any product.

## Note to Purchaser: Limited License

KAPA Hyper Prep Kits are developed, designed and sold exclusively for research purposes and *in vitro* use. Neither the product, nor any individual component, has been tested for use in diagnostics or for drug development, nor is it suitable for administration to humans or animals. Please refer to the MSDS, which is available on request.

Certain applications of this product are covered by patents issued to parties other than Kapa Biosystems and applicable in certain countries. Purchase of this product does not include a license to perform any such applications. Users of this product may therefore be required to obtain a patent license depending upon the particular application and country in which the product is used.

Licensed under U.S. Patent nos. 5,338,671 and 5,587,287 and corresponding patents in other countries.

Agencourt, AMPure and SPRI are registered trademarks of Beckman Coulter, Inc. Illumina and TruSeq are trademarks or registered trademarks of Illumina, Inc. xGen and Lockdown are trademarks of Integrated DNA Technologies. Qubit and PicoGreen are registered trademarks of Life Technologies. LabChip is a registered trademark of PerkinElmer, Inc. Nimblegen is a trademark of Roche.