# KAPA Hyper Prep Kit 取扱説明書

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## KAPA Hyper Prep Kit

KR0961 – v5.16

This Technical Data Sheet provides product information and a detailed protocol for the KAPA Hyper Prep Kit.

### Contents

Product Description	2
Product Applications	2
Product Specifications	2
Shipping and Storage	2
Handling	2
Quality Control	2
<b>Important Parameters</b>	<b>3</b>
Automated Library Construction	3
Safe Stopping Points	3
Input DNA and Fragmentation	3
Adapter Design and Concentration	4
Post-ligation Processing	5
Reaction Cleanups	5
Size Selection	6
Library Amplification	6
Evaluating the Success of Library Construction	8
<b>Process Workflow</b>	<b>10</b>
<b>Library Construction Protocol</b>	<b>11</b>
Appendix 1: Size Selection	13
Appendix 2: Library Construction Guidelines for the Roche® NimbleGen™ SeqCap™ EZ Target Capture System	15
Restrictions and Liabilities	16
Note to Purchaser: Limited Product Warranty	16
Note to Purchaser: Limited License	16

Kapa/Roche Kit Codes and Components		
<b>KK8500</b> 07962312001	End Repair & A-Tailing Buffer	60 µL
	End Repair & A-Tailing Enzyme	25 µL
<b>KK8501*</b> 07962339001	Ligation Buffer	245 µL
	DNA Ligase	80 µL
	KAPA HiFi HotStart ReadyMix (2X)*	250 µL
8 libraries	Library Amplification Primer Mix (10X)*	50 µL
<b>KK8502</b> 07962347001	End Repair & A-Tailing Buffer	210 µL
	End Repair & A-Tailing Enzyme	90 µL
<b>KK8503*</b> 07962355001	Ligation Buffer	900 µL
	DNA Ligase	300 µL
	KAPA HiFi HotStart ReadyMix (2X)*	690 µL
24 libraries	Library Amplification Primer Mix (10X)*	138 µL
<b>KK8504</b> 07962363001	End Repair & A-Tailing Buffer	930 µL
	End Repair & A-Tailing Enzyme	400 µL
<b>KK8505*</b> 07962371001	Ligation Buffer	3.8 mL
	DNA Ligase	1.26 mL
	KAPA HiFi HotStart ReadyMix (2X)*	3.0 mL
96 libraries	Library Amplification Primer Mix (10X)*	600 µL

\*KK8501, KK8503, and KK8505 are available for PCR-free workflows, and do not contain library amplification reagents.

### Quick Notes

- This kit provides a versatile, streamlined DNA library construction protocol. Libraries for Illumina® sequencing may be prepared from a wide range of DNA samples and inputs (1 ng – 1 µg) in 2 – 3 hrs, with minimal hands-on time. DNA input is appropriately fragmented double-stranded DNA.
- The novel one-tube chemistry improves library yield and quality, particularly for FFPE and low-input samples.
- The protocol is easy to automate. Generous reagent excesses are supplied in 96-reaction kits to accommodate the dead volumes required for automated liquid handling.
- Kits contain all the reagents for library construction and high-efficiency, low-bias library amplification, except for adapters and beads. KAPA Pure Beads (KK8000, KK8001, KK8002) and KAPA Adapters are sold separately. Kits without an amplification module are available for PCR-free workflows.
- The **Process Workflow** provides an overview of the library construction process. **Appendix 1** provides protocols for bead-based, double-sided size selection. **Appendix 2** provides details for integration with the SeqCap EZ workflow.

## Product Description

The KAPA Hyper Prep Kit provides a versatile, streamlined protocol for the rapid construction of libraries for Illumina® sequencing from fragmented, double-stranded DNA (dsDNA). The novel chemistry and streamlined, one-tube protocol improves the efficiency and consistency of library construction across a wide range of sample types and inputs.

The workflow combines enzymatic steps and employs minimal bead-based cleanups, thereby reducing sample handling and overall library preparation time to 2 – 3 hrs. The kit contains all of the enzymes and reaction buffers required for:

1. end repair and A-tailing, which produces end-repaired, 5'-phosphorylated, 3'-dA-tailed dsDNA fragments;
2. adapter ligation, during which dsDNA adapters with 3'-dTTP overhangs are ligated to 3'-dA-tailed molecules;
3. library amplification (optional), which employs high-fidelity, low-bias PCR to amplify library fragments carrying appropriate adapter sequences on both ends.

The kit provides a single, concentrated buffer and a single enzyme mixture for each of the two library construction steps. This offers the best combination of product stability, convenience and efficiency. Adapters and beads required for cleanups after adapter ligation and library amplification are not included. KAPA Pure Beads (KK8000, KK8001, KK8002) and KAPA Adapters are sold separately.

In order to maximize sequence coverage uniformity, it is critical to minimize library amplification bias. KAPA HiFi DNA Polymerase is designed for low-bias, high-fidelity PCR, and is the reagent of choice for NGS library amplification.<sup>1,2,3,4</sup> KAPA Hyper Prep Kits include KAPA HiFi HotStart ReadyMix (2X), a ready-to-use PCR mix comprising all the components for library amplification—except primers and template. Kits also include Library Amplification Primer Mix (10X), designed for the high-efficiency amplification of Illumina libraries flanked by adapters containing the P5 and P7 flow cell sequences. Kits without the amplification module (KK8501, KK8503, and KK8505) are available for PCR-free workflows. They may also be combined with KAPA HiFi Real-time Library Amplification Kits (KK2701 and KK2702), or with KAPA HiFi HotStart Uracil+ ReadyMix (KK2801 and KK2802) for the amplification of bisulfite-converted libraries.

1. Oyola, S.O., et al., *BMC Genomics* **13**, 1 (2012).
2. Quail, M.A., et al., *Nature Methods* **9**, 10 (2012).
3. Quail, M.A., et al., *BMC Genomics* **13**, 341 (2012).
4. Ross, M.G., et al., *Genome Biology* **14**, R51 (2013).

## Product Applications

KAPA Hyper Prep Kits are ideally suited for low- and high-throughput NGS library construction workflows that require end repair, A-tailing, adapter ligation and library amplification (optional). Kits are designed for library construction from a wide range of sample types and

inputs (1 ng – 1 µg), and are compatible with complex genomic DNA, cell-free/circulating tumor DNA and low-quality DNA such as FFPE samples. For small genomes, cell-free/circulating tumor DNA and lower complexity samples such as ChIP DNA, amplicons or cDNA (for RNA-seq), library construction from lower inputs (~100 pg or more) may be attempted.

The protocol is automation-friendly and may be incorporated into workflows for a wide range of NGS applications, including:

- whole-genome, shotgun sequencing
- whole exome or targeted sequencing, using Roche® NimbleGen™ SeqCap™ EZ, Agilent SureSelect, Illumina TruSeq®, or IDT xGen™ Lockdown™ Probes, or other hybridization capture systems
- ChIP-seq
- RNA-seq (starting with cDNA)
- methyl-seq (in combination with KAPA HiFi HotStart Uracil+ ReadyMix for library amplification)

## Product Specifications

### Shipping and Storage

The enzymes provided in this kit are temperature sensitive, and appropriate care should be taken during shipping and storage. KAPA Hyper Prep Kits are shipped on dry ice or ice packs, depending on the destination country. Upon receipt, immediately store enzymes and reaction buffers at -15°C to -25°C in a constant-temperature freezer. When stored under these conditions and handled correctly, the kit components will retain full activity until the expiry date indicated on the kit label.

### Handling

Always ensure that KAPA Hyper Prep Kit components have been fully thawed and thoroughly mixed before use. The KAPA Hyper Prep End Repair & A-tailing Buffer and Ligation Buffer may contain precipitates when thawed at 4°C. These buffers must be thawed at room temperature and vortexed thoroughly before use. KAPA HiFi HotStart ReadyMix (2X) contains isostabilizers and may not freeze completely, even when stored at -15°C to -25°C. Nevertheless, always ensure that the ReadyMix is fully thawed and thoroughly mixed before use. Reaction master mixes prepared from the enzymes and buffers for end repair and A-tailing, as well as for ligation, are very viscous and require special attention during pipetting. Keep all enzyme components and master mixes on ice as long as possible during handling and preparation.

### Quality Control

All kit components are subjected to stringent functional quality control, are free of detectable contaminating exo- and endonuclease activities, and meet strict requirements with respect to DNA contamination. Please contact Technical Support at [kapabiosystems.com/support](http://kapabiosystems.com/support) for more information.

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
## Recommended QC Metrics

## Library Construction Protocol

**Note:** This protocol does not include size selection. Please refer to **Appendix 1** for detailed double-sided size selection protocols.

### 1. End Repair and A-tailing

- 1.1 Assemble each end repair and A-tailing reaction in a tube or well of a PCR plate as follows:

Component	Volume
Fragmented, double-stranded DNA	50 µL
End Repair & A-Tailing Buffer*	7 µL
End Repair & A-Tailing Enzyme Mix*	3 µL
<b>Total volume:</b>	<b>60 µL</b>

\*The buffer and enzyme mix should preferably be pre-mixed and added in a single pipetting step. Premixes are stable for ≤24 hrs at room temperature, for ≤3 days at 4°C, and for ≤4 weeks at -20°C.

- 1.2 Vortex gently and spin down briefly. Return the plate/tube(s) to ice. Proceed immediately to the next step.
- 1.3 Incubate in a thermocycler programmed as outlined below:

Step	Temp	Time
End repair and A-tailing	20°C	30 min
	65°C*	30 min
HOLD	4°C**	∞

\*A heated lid is required for this incubation. If possible, set the temperature of the lid at 85°C, instead of the usual ~105°C.

\*\*If proceeding to the adapter ligation reaction setup without any delay, the reaction may be cooled to 20°C instead of 4°C.

- 1.4 Proceed immediately to **Adapter Ligation** (step 2).

### 2. Adapter Ligation

- 2.1 Dilute adapter stocks to the appropriate concentration, as outlined in Table 2 (p. 5).
- 2.2 In the same plate/tube(s) in which end repair and A-tailing was performed, assemble each adapter ligation reaction as follows:

Component	Volume
End repair and A-tailing reaction product	60 µL
Adapter stock (concentration as required)	5 µL
PCR-grade water*	5 µL
Ligation Buffer*	30 µL
DNA Ligase*	10 µL
<b>Total volume:</b>	<b>110 µL</b>

\*The water, buffer and ligase enzyme should preferably be premixed and added in a single pipetting step. Premixes are stable for ≤24 hrs at room temperature, for ≤3 days at 4°C, and for ≤4 weeks at -20°C.

- 2.3 Mix thoroughly and centrifuge briefly.

- 2.4 Incubate at 20°C for 15 min.

**Note:** to achieve higher conversion rates and library yields, particularly for low-input samples, consider increasing the ligation time—to a maximum of 4 hrs at 20°C, or overnight at 4°C. Please note that longer ligation times may lead to increased levels of adapter-dimer. Adapter concentrations may have to be optimized if ligation times are extended significantly.

- 2.5 Proceed immediately to the next step.

### 3. Post-ligation Cleanup

- 3.1 In the same plate/tube(s), perform a 0.8X bead-based cleanup by combining the following:

Component	Volume
Adapter ligation reaction product	110 µL
KAPA Pure Beads	88 µL
<b>Total volume:</b>	<b>198 µL</b>

- 3.2 Mix thoroughly by vortexing and/or pipetting up and down multiple times.
- 3.3 Incubate the plate/tube(s) at room temperature for 5 – 15 min to bind DNA to the beads.
- 3.4 Place the plate/tube(s) on a magnet to capture the beads. Incubate until the liquid is clear.
- 3.5 Carefully remove and discard the supernatant.
- 3.6 Keeping the plate/tube(s) on the magnet, add 200 µL of 80% ethanol.
- 3.7 Incubate the plate/tube(s) on the magnet at room temperature for ≥30 sec.
- 3.8 Carefully remove and discard the ethanol.
- 3.9 Keeping the plate/tube(s) on the magnet, add 200 µL of 80% ethanol.
- 3.10 Incubate the plate/tube(s) on the magnet at room temperature for ≥30 sec.
- 3.11 Carefully remove and discard the ethanol. Try to remove all residual ethanol without disturbing the beads.
- 3.12 Dry the beads at room temperature for 3 – 5 min, or until all of the ethanol has evaporated. **Caution: over-drying the beads may result in reduced yield.**
- 3.13 Remove the plate/tube(s) from the magnet.
- 3.14 Thoroughly resuspend the beads:
- in 25 µL of elution buffer (10 mM Tris-HCl, pH 8.0 – 8.5) to proceed with **Library Amplification** (step 4), or
  - in 55 µL of elution buffer (10 mM Tris-HCl, pH 8.0 – 8.5) to proceed with double-sided size selection (**Appendix 1**).

- 3.15 Incubate the plate/tube(s) at room temperature for 2 min to elute DNA off the beads.
- 3.16 Place the plate/tube(s) on a magnet to capture the beads. Incubate until the liquid is clear.
- 3.17 Transfer the clear supernatant to a new plate/tube(s):
  - to proceed with **Library Amplification** (step 4), transfer 20  $\mu\text{L}$  of supernatant, or
  - to proceed with double-sided size selection (**Appendix 1**), transfer 50  $\mu\text{L}$  of supernatant.

## 4. Library Amplification

Note: Please refer to **Important Parameters: Library Amplification** (pp. 6 – 7) and the **KAPA NGS Library Preparation Technical Guide** for more information on optimizing library amplification.

- 4.1 Assemble each library amplification reaction as follows:

Component	Volume
KAPA HiFi HotStart ReadyMix (2X)	25 $\mu\text{L}$
KAPA Library Amplification Primer Mix (10X)*	5 $\mu\text{L}$
Adapter-ligated library	20 $\mu\text{L}$
<b>Total volume:</b>	<b>50 <math>\mu\text{L}</math></b>

\*Or another, suitable 10X library amplification primer mix. The recommended final concentration of each primer in the library amplification reaction is 0.5 – 4  $\mu\text{M}$ . Also refer to **Important Parameters: Library Amplification** (p. 6).

- 4.2 Mix thoroughly and centrifuge briefly.
- 4.3 Amplify using the following cycling protocol:

Step	Temp	Duration	Cycles
Initial denaturation	98°C	45 sec	1
Denaturation	98°C	15 sec	Minimum number required for optimal amplification (Table 3 or 4)
Annealing*	60°C	30 sec	
Extension	72°C	30 sec	
Final extension	72°C	1 min	1
HOLD	4°C	$\infty$	1

\*Optimization of the annealing temperature may be required for non-standard (i.e., other than Illumina® TruSeq®) adapter/primer combinations.

- 4.4 Proceed directly to **Post-amplification Cleanup** (step 5).

## 5. Post-amplification Cleanup

- 5.1 In the library amplification plate/tube(s) perform a 1X bead-based cleanup by combining the following:

Component	Volume
Library amplification reaction product	50 $\mu\text{L}$
KAPA Pure Beads	50 $\mu\text{L}$
<b>Total volume:</b>	<b>100 <math>\mu\text{L}</math></b>

- 5.2 Mix thoroughly by vortexing and/or pipetting up and down multiple times.
- 5.3 Incubate the plate/tube(s) at room temperature for 5 – 15 min to bind DNA to the beads.
- 5.4 Place the plate/tube(s) on a magnet to capture the beads. Incubate until the liquid is clear.
- 5.5 Carefully remove and discard the supernatant.
- 5.6 Keeping the plate/tube(s) on the magnet, add 200  $\mu\text{L}$  of 80% ethanol.
- 5.7 Incubate the plate/tube(s) on the magnet at room temperature for  $\geq 30$  sec.
- 5.8 Carefully remove and discard the ethanol.
- 5.9 Keeping the plate/tube(s) on the magnet, add 200  $\mu\text{L}$  of 80% ethanol.
- 5.10 Incubate the plate/tube(s) on the magnet at room temperature for  $\geq 30$  sec.
- 5.11 Carefully remove and discard the ethanol. Try to remove all residual ethanol without disturbing the beads.
- 5.12 Dry the beads at room temperature for 3 – 5 min, or until all of the ethanol has evaporated. **Caution: over-drying the beads may result in reduced yield.**
- 5.13 Remove the plate/tube(s) from the magnet.
- 5.14 Thoroughly resuspend the beads in an appropriate volume of elution buffer (10 mM Tris-HCl, pH 8.0 – 8.5) or PCR-grade water. **Always use PCR-grade water if proceeding to target capture.**  
**Note:** If proceeding with a second post-amplification cleanup, or double-sided size selection (**Appendix 1**), resuspend the beads in 55  $\mu\text{L}$  of elution buffer.
- 5.15 Incubate the plate/tube(s) at room temperature for 2 min to elute DNA off the beads.
- 5.16 Place the plate/tube(s) on a magnet to capture the beads. Incubate until the liquid is clear.
- 5.17 Transfer the clear supernatant to a new plate/tube(s) and proceed with size selection (refer to **Appendix 1**), library QC, target capture or sequencing, as appropriate. Store purified, amplified libraries at 4°C for 1 – 2 weeks, or at -20°C.

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## Appendix 2: Library Construction Guidelines for the Roche® NimbleGen™ SeqCap™ EZ Target Capture System

The KAPA Hyper Prep Kit may be used (instead of a KAPA HTP/LTP Library Preparation Kit) for the construction of pre-capture libraries in the SeqCap EZ workflow—to ensure a more streamlined protocol, shorter turnaround times and improved performance with low-input and challenging samples, such as FFPE and cell-free/free circulating tumor DNA.

### Library Construction

Detailed instructions for specific steps or processes are included below.

#### DNA Input and Fragmentation

- The standard SeqCap EZ Library SR workflow specifies 100 ng of Covaris-sheared DNA as the input into library construction (the End Repair and A-tailing reaction).

**Note:** libraries may be constructed from lower inputs, but sequencing metrics are not guaranteed. To use the KAPA Hyper Prep Kit in a SeqCap EZ workflow with lower DNA inputs, please consult **Important Parameters** (pp. 3 – 9) or contact Technical Support at [kapabiosystems.com/support](http://kapabiosystems.com/support) for more information.

- Shear the input DNA to an average size of 180 – 220 bp using a Covaris instrument. If an optimized protocol is not available then start with the guidelines provided by the supplier of the Covaris instrument and fine-tune the parameters to yield the appropriate DNA fragmentation size distribution.

#### End Repair and A-tailing

- Follow instructions outlined in the **Library Construction Protocol: End Repair and A-tailing** (step 1).

#### Adapter Ligation

- Follow instructions outlined in the **Library Construction Protocol: Adapter Ligation** (step 2).
- Use 5 µL of undiluted SeqCap Indexed Adapter per reaction.

#### Post-ligation Processing

- Perform the 0.8X post-ligation cleanup as outlined in the **Library Construction Protocol: Post-ligation Cleanup** (step 3). Resuspend the beads with adapter-ligated DNA in 55 µL elution buffer, and recover 50 µL of supernatant.
- Perform a 0.7X – 0.9X double-sided size selection, with KAPA Pure Beads as outlined in **Appendix 1** (p. 13).
- Resuspend beads with size-selected, adapter-ligated DNA in 22 – 25 µL elution buffer, and recover 20 µL of supernatant.

#### Pre-capture Library Amplification (LM-PCR)

- Amplify 20 µL of each pre-capture library, using the reaction setup and cycling parameters outlined in **Library Construction Protocol: Library Amplification** (step 4).
- Use 5 µL of the KAPA Library Amplification Primer Mix (10X) supplied in the KAPA Hyper Prep Kit for each amplification reaction.
- Perform 9 cycles of pre-capture amplification.

**Note:** 9 cycles of pre-capture amplification should yield  $\geq 1$  µg of pre-capture library (i.e., a sufficient amount to perform individual captures). If you are pooling libraries before capture, the amount of DNA needed per capture is  $\geq 1$  µg/n, where n = the number of libraries per capture. The number of pre-capture amplification cycles may therefore be reduced to achieve the lowest possible duplication rates and highest library complexity prior to capture.

#### Cleanup after Pre-capture Library Amplification (LM-PCR)

- Perform a 1X bead-based cleanup of the pre-capture LM-PCR products, as outlined in **Library Construction Protocol: Post-amplification Cleanup** (step 5).

#### Quality Control of Pre-capture Libraries

- Please refer to the SeqCap EZ Library SR User's Guide v5.1 (Chapter 4, Step 5) (Roche document number 06588786001, 09/15).

#### Target Enrichment

- Please refer to Chapters 5 – 8 of the SeqCap EZ Library SR User's Guide v5.1 to complete the workflow and generate sequencing-ready, enriched libraries.

## Restrictions and Liabilities

This technical data sheet is provided “as is” and Kapa Biosystems assumes no responsibility for any typographical, technical, or other inaccuracies. The document is subject to change, without notice, in future editions.

To the maximum extent permitted by applicable law, Kapa Biosystems disclaims all warranties, either express or implied, with regard to this technical data sheet and any information contained herein, including but not limited to the implied warranties of merchantability and fitness for a particular purpose. Kapa Biosystems shall not be liable for errors or for incidental or consequential damages in connection with the furnishing, use, or performance of this document or of any information contained herein.

This document might contain references to third party sources of information, hardware or software, products, or services and/or third party web sites (collectively the “Third-Party Information”). Kapa Biosystems does not control, and is not responsible for, any Third-Party Information. The inclusion of Third-Party Information in this document does not imply endorsement by Kapa Biosystems of the Third-Party Information or the third party in any way.

Kapa Biosystems is not responsible nor will be liable in any way for your use of any software or equipment that is not supplied by Kapa Biosystems in connection with your use of Kapa Biosystems products.

Kapa Biosystems does not in any way guarantee or represent that you will obtain satisfactory results from using Kapa Biosystems products as described herein. The only warranties provided to you are included in the Limited Warranty enclosed with this document. You assume all risk in connection with your use of Kapa Biosystems products.

## Note to Purchaser: Limited Product Warranty

Any product that does not meet the performance standards stated in the product specification sheet will be replaced at no charge. This warranty limits our liability to the replacement of the product. No other warranties of any kind, express or implied, including without limitation, implied warranties of merchantability or fitness for a particular purpose, are provided by Kapa Biosystems. Kapa Biosystems shall have no liability for any direct, indirect, consequential or incidental damages arising out of the use, the results of use or the inability to use any product.

## Note to Purchaser: Limited License

KAPA Hyper Prep Kits are developed, designed and sold exclusively for research purposes and *in vitro* use. Neither the product, nor any individual component, has been tested for use in diagnostics or for drug development, nor is it suitable for administration to humans or animals. Please refer to the MSDS, which is available on request.

Certain applications of this product are covered by patents issued to parties other than Kapa Biosystems and applicable in certain countries. Purchase of this product does not include a license to perform any such applications. Users of this product may therefore be required to obtain a patent license depending upon the particular application and country in which the product is used.

Licensed under U.S. Patent nos. 5,338,671 and 5,587,287 and corresponding patents in other countries.