



Desalting and Buffer Exchange with Vivaspin® Centrifugal Concentrators



#04

Application
Note

#05

#06

#07

#08

Pieter Eyckermann, Rik McRae,
Andreas Kocourek and Robert Zeidler

Introduction

Vivaspin® centrifugal concentrators, with patented vertical membrane technology, combine fast filtration with high recovery of target proteins. This makes Vivaspin® the technology of choice for desalting or buffer exchange, avoiding lengthy dialysis steps.

While proteins are retained by an appropriate ultrafiltration membrane, salts can pass freely through, independent of protein concentration or membrane MWCO. In consequence, the composition of the buffer in the flow-through and retentate is unchanged after protein concentration. By diluting the concentrate back to the original volume, the salt concentration is lowered. The concentrate can be diluted with water or salt-free buffer if simple desalting is required; however, it is also possible to dilute the concentrate with a new buffer, thereby exchanging the buffering substance entirely. For example, a 10 ml protein sample containing 500 mM salt, if concentrated 100x still contains 500 mM salt. If this concentrate is then diluted 100x with water or salt-free buffer, the protein concentration returns to normal, while the salt concentration is reduced 100x to only 5 mM, (i.e. a 99% reduction in salt).

The protein sample can then be concentrated again to the desired level, or the buffer exchange can be repeated to reduce the salt concentration even further before a final concentration of the protein. This process is called 'diafiltration'. For proteins with a tendency to precipitate at higher concentrations, it is possible to perform several diafil-

tration steps in sequence, with the protein concentrated each time to only 5 or 10x. For example, if a precipitous protein sample is concentrated to 5x then diluted back to the original volume, and this process is repeated a further two times, this still results in a >99% reduction in salt concentration, without over concentrating the protein.

Desalting and Buffer Exchange Procedure (See Figure 1.)

1. Select the most appropriate MWCO for your sample. For maximum recovery, select a MWCO $\frac{1}{2}$ to $\frac{1}{3}$ the molecular size of the species of interest.
2. Fill concentrator with up to the maximum volume stated in the device operating instructions*, (e.g. 20 ml if Vivaspin® 20 is used).
3. If the sample is smaller than the maximum device volume*, it can be diluted up to the maximum volume before the first centrifugation step. This will help increase the salt removal rate.
4. Centrifuge for the recommended amount of time at an appropriate spin speed for your Vivaspin® model*.
5. Empty filtrate container†.
6. Refill concentrator with an appropriate solvent.
7. Centrifuge again as before.
8. Empty filtrate container†.
9. Recover the concentrated, de-salted sample from the bottom of the concentrate pocket with a pipette.

Notes

* For guidance on maximum fill volumes, spin speeds and suggested spin times, please refer to the Operating Instructions that accompany your Vivaspin® products.

† Filtrate volumes should be retained until the concentrated sample has been analyzed.




Figure 1: Step-by-step method for desalting and concentration

Test Results

As the results below show, the efficient design of Vivaspin® devices allowed >95% of the salt to be removed during the first centrifugation step. Only one subsequent centrifugation step was needed to increase the typical salt removal to 99% with >92% recovery of the sample.

Vivaspin® 20

	MWCO 5 kDa		30 kDa		50 kDa		100 kDa	
	Protein Recovery	NaCl Removal	Protein Recovery	NaCl Removal	Protein Recovery	NaCl Removal	Protein Recovery	NaCl Removal
	Cytochrome C 0.25 mg/ml		BSA 1 mg/ml		BSA 1 mg/ml		IgG 1 mg/ml	
Spin 1	100 %	99 %	97 %	99 %	97 %	99 %	90 %	98 %
Spin 2	96 %	100 %	92 %	100 %	93 %	100 %	87 %	100 %

Four Vivaspin® 20 devices of each cut-off were tested with 20 ml of solution. Each of the solutions contained 500 mM NaCl. Each spin was performed at 4,000 × g. The devices > 5kDa were spun for 30 min. The devices with 5 kDa were spun 45 min. After the first and second spin, the retentate was brought up to 20 ml with ultra pure water from the Arium system (Sartorius). OD readings were taken at 410 nm for the Cytochrome C and 280 nm for the BSA and IgG samples. Salt concentration was measured with a Qcond 2200 conductivity measuring instrument.

Vivaspin® 6

	MWCO 5 kDa		30 kDa		50 kDa		100 kDa	
	Protein Recovery	NaCl Removal	Protein Recovery	NaCl Removal	Protein Recovery	NaCl Removal	Protein Recovery	NaCl Removal
	Cytochrome C 0.25 mg/ml		BSA 1 mg/ml		BSA 1 mg/ml		IgG 1 mg/ml	
Spin 1	98 %	99 %	92 %	99 %	93 %	99 %	92 %	98 %
Spin 2	85 %	100 %	86 %	100 %	83 %	100 %	89 %	100 %

Four Vivaspin® 6 devices of each cut-off were tested with 6 ml of solution. Each of the solutions contained 500 mM NaCl. Each spin was performed at 4,000 × g. The devices > 5 kDa were spun for 30 min. The devices with 5 kDa were spun 45 min. After the first and the second spin the retentate was brought up to 6 ml with ultra pure water from the Arium system (Sartorius) OD readings were taken at 410 nm for the Cytochrome C and 280 nm for the BSA and IgG samples. Salt concentration was measured with a Qcond 2200 conductivity measuring instrument.

Sales and Service Contacts

For further contacts, visit www.sartorius-stedim.com

Europe

Germany

Sartorius Stedim Biotech GmbH
August-Spindler-Strasse 11
37079 Goettingen

Phone +49.551.308.0
Fax +49.551.308.3289

www.sartorius-stedim.com

Sartorius Stedim Systems GmbH
Robert-Bosch-Strasse 5-7
34302 Guxhagen

Phone +49.5665.407.0
Fax +49.5665.407.2200

France

Sartorius Stedim Biotech S.A.
ZI Les Paluds
Avenue de Jouques - BP 1051
13781 Aubagne Cedex

Phone +33.442.845600
Fax +33.442.845619

Sartorius Stedim France SAS
ZI Les Paluds
Avenue de Jouques - CS 71058
13781 Aubagne Cedex

Phone +33.442.845600
Fax +33.442.846545

Austria

Sartorius Stedim Austria GmbH
Franzosengraben 12
1030 Vienna

Phone +43.1.7965763.18
Fax +43.1.796576344

Belgium

Sartorius Stedim Belgium N.V.
Leuvensesteenweg, 248/B
1800 Vilvoorde

Phone +32.2.756.06.80
Fax +32.2.756.06.81

Denmark

Sartorius Stedim Nordic A/S
Hoerskaetten 6D, 1.
2630 Taastrup

Phone +45.7023.4400
Fax +45.4630.4030

Hungary

Sartorius Stedim Hungária Kft
Kagyló u. 5
2092 Budakeszi

Phone +36.23.457.227
Fax +36.23.457.147

Italy

Sartorius Stedim Italy S.p.A.
Via dell'Antella, 76/A
50012 Antella-Bagno a Ripoli (FI)

Phone +39.055.63.40.41
Fax +39.055.63.40.526

Netherlands

Sartorius Stedim Netherlands B.V.
Edisonbaan 24
3439 MN Nieuwegein

Phone +31.30.6025080
Fax +31.30.6025099

Poland

Sartorius Stedim Poland Sp. z o.o.
ul. Wrzesinska 70
62-025 Kostrzyn

Phone +48.61.647.38.40
Fax +48.61.879.25.04

Russian Federation

000 „Sartorius ICR“ and 000 „Biohit“
Uralskaya str. 4, Lit. B
199155, Saint-Petersburg

Phone +7.812.327.5.327
Fax +7.812.327.5.323

Spain

Sartorius Stedim Spain SA
C/Isabel Colbrand 10,
Oficina 70
Poligono Industrial de Fuencarral
28050 Madrid

Phone +34.90.2110935
Fax +34.91.3589623

Switzerland

Sartorius Stedim Switzerland AG
Ringstr. 24 a
8317 Tagelswangen

Phone +41.52.354.36.36
Fax +41.52.354.36.46

U.K.

Sartorius Stedim UK Limited
Longmead Business Park
Blenheim Road, Epsom
Surrey KT19 9 QQ

Phone +44.1372.737159
Fax +44.1372.726171

America

USA

Sartorius Stedim North America Inc.
5 Orville Drive, Suite 200
Bohemia, NY 11716

Toll-Free +1.800.368.7178
Fax +1.631.254.4253

Argentina

Sartorius Argentina S.A.
Int. A. Avalos 4251
B1605ECS Munro
Buenos Aires

Phone +54.11.4721.0505
Fax +54.11.4762.2333

Brazil

Sartorius do Brasil Ltda
Av. Dom Pedro I, 241
Bairro Vila Pires
Santo André
São Paulo
Cep 09110-001

Phone +55.11.4451.6226
Fax +55.11.4451.4369

Mexico

Sartorius de México S.A. de C.V.
Circuito Circunvalación Poniente No. 149
Ciudad Satélite
53100, Estado de México
México

Phone +52.5555.62.1102
Fax +52.5555.62.2942

Asia | Pacific

Australia

Sartorius Stedim Australia Pty. Ltd.
Unit 5, 7-11 Rodeo Drive
Dandenong South Vic 3175

Phone +61.3.8762.1800
Fax +61.3.8762.1828

China

Sartorius Stedim Biotech (Beijing)
Co. Limited
Airport Industrial Zone B
No. 33 Yu'an Road
Beijing 101300, Shunyi District

Phone +86.10.80426516
Fax +86.10.80426580

Sartorius Stedim Biotech (Beijing)
Co. Limited
Shanghai Branch office
Room 618, Tower 1, German Centre,
Shanghai, PRC., 201203

Phone +86.21.28986393
Fax +86.21.28986392.11

Sartorius Stedim Biotech (Beijing)
Co. Limited
Guangzhou representative office
Room 704, Broadway Plaza,
No. 233-234 Dong Feng West Road
Guangzhou 510180

Phone +86.20.8351.7921
Fax +86.20.8351.7931

India

Sartorius Stedim India Pvt. Ltd.
#69/2-69/3, Jakkasandra
Kunigal Road, Nelamangala Tq
Bangalore - 562 123

Phone +91.80.4350.5248
Fax +91.80.4350.5253

Japan

Sartorius Stedim Japan K.K.
Kiba Park Bldg
5-11-13 Kiba
Koto-ku
Tokyo 135-0042

Phone +81.3.5639.9981
Fax +81.3.5639.9983

Malaysia

Sartorius Stedim Malaysia Sdn. Bhd.
Lot L3-E-3B, Enterprise 4
Technology Park Malaysia
Bukit Jalil
57000 Kuala Lumpur, Malaysia

Phone +60.3.8996.0622
Fax +60.3.8996.0755

Singapore

Sartorius Stedim Singapore Pte. Ltd.
1 Science Park Road,
The Capricorn, #05-08A,
Singapore Science Park 2
Singapore 117528

Phone +65.6872.3966
Fax +65.6778.2494

South Korea

Sartorius Korea Biotech Co., Ltd.
8th Floor, Solid Space B/D,
PanGyoYeok-Ro 220, BunDang-Gu
SeongNam-Si, GyeongGi-Do, 463-400

Phone +82.31.622.5700
Fax +82.31.622.5799



◀▶ www.sartorius-stedim.com