

# MINIGRIP

A Member of Inteplast Group

# VWR

1.800.932.5000

vwr.com

# Lab Guard®

## DESTROYABLE

### Biohazard Symbol



Open bag at patented TearZone® area -- The biohazard symbol is destroyed!



Throw the uncontaminated bag remains away in your **REGULAR WASTE!** If the bag remains are contaminated, throw them away in **RED** waste.\*



FDA Registered.



Product Description	Case Size / Packaging	Minigrip Part #	VWR Cat. No.
6 x 9, 3-Wall, TearZone® Feature	1000 (20 packs of 50 ea) DispenserBag® Pkg	SBL2X69B	56766-370
6 x 9, 3-Wall, TearZone® Feature	1000 (2 - 500 pack cartons) Dispenser Box Pkg	SBL2X69B5	95059-670



VWR 1.800.932.5000  
vwr.com

Prices and product details are current when published; subject to change without notice. | Certain products may be limited by federal, state, provincial, or local regulations. | VWR makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC. All prices are in US dollars unless otherwise noted. Offers valid in US and Canada, void where prohibited by law or company policy, while supplies last. | VWR, the VWR logo and variations on the foregoing are registered (®) or unregistered trademarks and service marks, of VWR International, LLC and its related companies. All other marks referenced are registered by their respective owner(s). | Visit vwr.com to view our privacy policy, trademark owners and additional disclaimers. ©2013 VWR International, LLC. All rights reserved.

\* Minigrip is not responsible for determining waste stream choice or if bag is contaminated.