



## Certificate of Quality

### Excella® Microcentrifuge Tubes

Article No.	1.5 mL	2.0 mL
Microcentrifuge Tube	9104021000	9104021001

Eppendorf certifies that this product is produced in a fully controlled environment that ensures the product is free from detectable RNase, DNase, DNA and Endotoxins. Selected samples of product from this production environment have been tested for the following parameters on a recurring basis.

- **RNase testing method:** RNA Digestion  
Limits:  $< 1 \times 10^{-9}$  Kunitz Units
- **DNase testing method:** DNA Digestion  
Limits:  $< 1 \times 10^{-7}$  Kunitz Units
- **Human DNA testing method:** Polymerase Chain Reaction  
Limits:  $< 30 \times 10^{-9}$  mg
- **Endotoxin testing method:** Kinetic-chromogenic LAL-Test  
limits:  $< 0.06$  EU/mL

Hamburg, March 2024

Mirko Blohm  
Head of Global Supplier Quality Assurance

Katlin Heinemeyer  
Business Development Manager

Eppendorf SE, Barkhausenweg 1, 22339 Hamburg, Germany  
Phone: + 49 40 53 801-0

[www.excella.bio](http://www.excella.bio)