



## Automation methods availability, by vendor

| Vendor          | Model   | KAPA DNA Library Preparation |                    |                                 |  | Roche Target Enrichment |                   | KAPA RNA Library Preparation |                        |                                       |                         |
|-----------------|---|------------------------------|--------------------|---------------------------------|--|-------------------------|-------------------|------------------------------|------------------------|---------------------------------------|-------------------------|
|                 |   | KAPA HyperPrep Kit           | KAPA HyperPlus Kit | KAPA Library Quantification Kit | KAPA Human Genomic DNA Quantification and QC Kit | Library Normalization   | HyperCap Workflow | SeqCap EZ Workflow           | KAPA RNA HyperPrep Kit | KAPA RNA HyperPrep Kit with RiboErase | KAPA mRNA HyperPrep Kit |
| PerkinElmer     | Sciclone G3 NGSx                                    | PerkinElmer                  | PerkinElmer        | PerkinElmer                     | PerkinElmer                                      |                         | PerkinElmer       | PerkinElmer                  | PerkinElmer            | PerkinElmer                           | PerkinElmer             |
|                 | Zephyr G3 NGS                                       | PerkinElmer                  | PerkinElmer        | PerkinElmer                     | PerkinElmer                                      |                         | PerkinElmer*      |                              | PerkinElmer            | PerkinElmer                           | PerkinElmer             |
|                 | JANUS G3R NGS Express                               | PerkinElmer                  | PerkinElmer        | PerkinElmer                     | PerkinElmer                                      |                         |                   |                              | PerkinElmer            | PerkinElmer                           |                         |
| Beckman Coulter | Biomek FX <sup>o</sup> and FX, dual pod             | Roche                        | Roche              | Roche                           | Roche*   |                         | Roche             |                              | Roche*                 | Roche*                                | Roche*                  |
|                 | Biomek NX <sup>o</sup> and NX <sup>o</sup> (Span-8) |                              |                    | Beckman Coulter                 |  |                         |                   |                              |                        |                                       |                         |
|                 | Biomek 4000   |                              |                    | Roche*                          |  |                         |                   |                              |                        |                                       |                         |
|                 | Biomek i5 Automated Workstation                     | Beckman Coulter*             | Beckman Coulter*   |                                 |  |                         |                   |                              |                        |                                       |                         |
| Agilent         | Biomek i7 Automated Workstation                     | Roche                        | Roche              |                                 |  |                         | Roche*            |                              |                        |                                       |                         |
|                 | Bravo NGS Workstation A                             | Agilent                      | Agilent            | Roche*                          |  |                         |                   | Agilent                      |                        |                                       |                         |
|                 | Bravo NGS Workstation B                             | Roche                        | Roche              | Roche                           | Roche*   | Roche*                  | Roche             | Agilent                      |                        |                                       | Agilent                 |
| Eppendorf       | epMotion 5075 TMX                                   | Eppendorf                    | Eppendorf          | Eppendorf                       |  |                         |                   |                              | Eppendorf*             |                                       | Eppendorf*              |
| Hamilton        | NGS STAR (8 channel)                                | Hamilton                     | Hamilton           | Hamilton                        |  |                         | Hamilton*         |                              | Hamilton*              | Hamilton*                             | Hamilton*               |
| Tecan           | Freedom Evo NGS Workstation                         | Tecan                        | Tecan              | Tecan*                          |  |                         |                   | Tecan                        | Tecan*                 | Tecan*                                | Tecan*                  |
| PrimaDiag       | ACSIA NGS (capture)                                 |                              | PrimaDiag          |                                 |  |                         |                   | PrimaDiag                    |                        |                                       |                         |
| Brand           | Liquid Handling System                              |                              |                    | Brand*                          |  |                         |                   |                              |                        |                                       |                         |

\*in development

Regional differences may apply, please contact your local Roche representative.

Published by:

**Roche Sequencing Solutions, Inc.**

4300 Hacienda Drive  
Pleasanton, CA 94588

[sequencing.roche.com](http://sequencing.roche.com)

For Research Use Only. Not for use in diagnostic procedures.  
HYPERCAP, KAPA and SEQCAP are trademarks of Roche. All other product names and trademarks are the property of their respective owners.  
©2018 Roche Sequencing Solutions, Inc. All rights reserved.