

Cat #: WPT10010 / WPT10011 / WPT10012 / WPT10200 / WPT10300 / WPT11000 / WPT11250
 WFT10010/WFT10011/WFT10012 / WFT10020 / WFT10100 / WFT10200 / WFT11000 / WFT11250

General Information

Description

Racked, Sterile, Universal Filtered Pipette Tips

Raw Material

100% Origin Polypropylene, USP Class VI, free of Heavy metal

Sterilization

E-Beam, SAL 10⁻⁶

Color

Natural

Compatibility

See Compatibility Chart

Features

Temperature

Approx. -80°C to +121°C

Autoclavability

Autoclave at 121°C/15psi for 20 mins

Shelflife

3 years after sterilization

Others

Free from Dnase, Rnase, Human DNA, Pyrogen, ATP, Endoxion

Packaging/Labeling

Packaging

96/rack 10 racks/ box 5 boxes/ carton

Lot Number

YYMMDD- Year Month Day

Dimension/Drawings

