

Spreader (#90050)

1. Description

Disposable spreader is used for spreading microorganisms, which is sterilized before the release.

2. Features

- ❖ For the spreading of microbe samples
- ❖ Autoclavable

3. Functional Test

Sterilized	:	Passed
RNase/DNase free	:	Passed
Human DNA free	:	Passed

4. Product Range

Cat. No.	Material	Sterile	Packaging
90050	PP	+	10 / 500

5. Certification


SPL Life Sciences hereby certifies that the product identified above, inspected to be in compliance with product quality specification and requirements as documented in our ISO 9001:2015 Quality Management System (K-QA-Q031478) in Korea.

6. Specifications & Technical Drawing

- Specifications

Dimensions	Material (Resin)	Packaging
Length _{grip} : 143.00 mm Length _{spreader} : 37.00 mm	Polypropylene (PP)	<ul style="list-style-type: none"> ▪ Pieces/Sleeve: 10 ▪ Pieces/Case: 500

- Technical Drawing



www.spllifesciences.com

For technical assistance, contact SPL R&D Center at:
 Tel: +82-31-533-4800; Fax: +82-31-533-1430; e-mail: spl@ispl.co.kr
 To place an order, contact your local distributor or
 Tel: +82-31-533-4800; Fax: +82-31-533-1430; e-mail: business@ispl.co.kr

SPL-TDS-SPR

For Research Use Only. Not intended for use in diagnostic or therapeutic procedures.
 © SPL Life Sciences Co., Ltd., Version 3, November 22th, 2018