

Test Tube

(#40005/40105/42105/41005/40205/40014/40014/40114/42114/41014)

1. Description

SPL Life Sciences provide a wide range of Test Tubes to meet various demands. Tubes are differentiated by volume, material or sterilization method.










2. Features

- ❖ Volumes: 5 or 14 ml
- ❖ Internal graduations
- ❖ Dual-position snap cap for easy handling
- ❖ No cap, non-sterile, for flow cytometry (Cat. No. 40205)
- ❖ Sterilized by E.O. gas (Cat. No. 40005, 40014, 41005, 41014)
- ❖ Sterilized by irradiation (Cat. No. 40105, 42105, 40114, 42114)
- ❖ Individually wrapped (Cat. No. 42105, 42114)
- ❖ Translucent polypropylene tubes with high chemical and thermal stability (Cat. No. 41005, 41014)

3. Functional Test

Sterilized	:	Passed
Non-cytotoxic	:	Passed
Non-pyrogenic	:	Passed
RNase / DNase - free	:	Passed
Human DNA - free	:	Passed

4. Product Range

TYPE	Cat. No.	Material (Tube/Cap)	Total Volume (ml)	External Dimension d x h (mm)	RCF rating	Sterile	Packaging
	40005	PS / LDPE	5.00	12.00 x 75.00	1400 x g	E.O.	25 / 500
	40105	PS / LDPE	5.00	12.00 x 75.00	1400 x g	Irradiation	25 / 500
	42105	PS / LDPE	5.00	12.00 x 75.00	1400 x g	Irradiation	1 / 500
	41005	PP / LDPE	5.00	12.00 x 75.00	3000 x g	E.O.	25 / 500
	40205	PS	5.00	12.00 x 75.00	1400 x g	-	100 / 1000
	40014	PS / LDPE	14.00	17.00 x 95.40	1400 x g	E.O.	25 / 500
	40114	PS / LDPE	14.00	17.00 x 95.40	1400 x g	Irradiation	25 / 500
	42114	PS / LDPE	14.00	17.00 x 95.40	1400 x g	Irradiation	1 / 500
	41014	PP / LDPE	14.00	17.00 x 95.40	3000 x g	E.O.	25 / 500

5. Certification

SPL Life Sciences hereby certifies that the product identified above, inspected to be in compliance with product quality specification and requirements as documented in our ISO 9001:2015 Quality Management System (K-QA-Q031478) in Korea.

SPL-TDS-TETU

For Research Use Only. Not intended for use in diagnostic or therapeutic procedures.
© SPL Life Sciences Co., Ltd., Version 3, May 12th, 2022

6. Product Photographs




7. Specifications & Technical Drawing

- Specifications


Dimensions		Material (Resin)		Packaging		
Cat.No. 40005/40105/42105/ 40205/41005	Cat.No. 40014/40114/ 42114/41014	Cat.No. 40005/40014/ 40105/42105/ 40114/42114/ 40205	Cat. No 41005/41014	Cat. No. 40205	Ca. No. 40005/40014/ 40105/40114/ 41005/41014	Cat. No. 42105/42114
Diameter: 17.70 mm Height: 76.60 mm	Diameter: 24.07 mm Height: 95.40 mm	Polystyrene (PS)	Polypropylene (PP)	▪ Pieces/Sleeve: 100 ▪ Pieces/Case: 1000	▪ Pieces/Sleeve: 25 ▪ Pieces/Case: 500	▪ Pieces/Sleeve: 1 ▪ Pieces/Case: 500

- Technical Drawing

Cat. No. 40005/40105/42105/40205/41005



Cat. No. 40014/40114/42114/41014



www.spllifesciences.com

For technical assistance, contact SPL R&D Center at:
Tel: +82-31-533-4800; Fax: +82-31-533-1430; e-mail: spl@ispl.co.kr
To place an order, contact your local distributor or
Tel: +82-31-533-4800; Fax: +82-31-533-1430; e-mail: business@ispl.co.kr

SPL-TDS-TETU

For Research Use Only. Not intended for use in diagnostic or therapeutic procedures.
© SPL Life Sciences Co., Ltd., Version 3, May 12th, 2022