

Conical Tube Rack

(#52015/52050)

1. Description

SPL Conical Tube Racks are economical and alternative solutions for stainless wire racks or other plastic racks. Racks are uniquely designed for handling and storage of 25 Conical Tubes. (15ml or 50ml respectively)

2. Features

- ❖ Excellent stackability
- ❖ Numeric labeling
- ❖ 3 colors: Blue, Green, Red

3. Product Range

Cat. No.	Material	Description	Packaging
52015	PP	25 Holes for 15 ml tubes	20
52050	PP	25 Holes for 50 ml tubes	20

4. Certification

SPL Life Sciences hereby certifies that the product identified above, inspected to be in compliance with product quality specification and requirements as documented in our ISO 9001:2015 Quality Management System (K-QA-Q031478) in Korea.

5. Product Photographs



SPL-TDS-COTURA

For Research Use Only. Not intended for use in diagnostic or therapeutic procedures.
© SPL Life Sciences Co., Ltd., Version 2, Sep19th, 2016


6. Specifications & Technical Drawing

- Specifications


Dimensions			Material (Resin)	Packaging	
52015	52050	52015		52050	
Length:	104.1 mm	Length:	Polypropylene (PP)	▪ Pieces/Case: 20	▪ Pieces/Case: 20
Width:	127.9 mm	Width:			
Height:	48.1 mm	Height:			

- Technical Drawing

Cat. No. 52015



Cat. No. 52050



www.spllifesciences.com

For technical assistance, contact SPL R&D Center at:

Tel: +82-31-533-4800; Fax: +82-31-533-1430; e-mail: spl@ispl.co.kr

To place an order, contact your local distributor or

Tel: +82-31-533-4800; Fax: +82-31-533-1430; e-mail: business@ispl.co.kr

SPL-TDS-COTURA

For Research Use Only. Not intended for use in diagnostic or therapeutic procedures.
© SPL Life Sciences Co., Ltd., Version 2, Sep19th, 2016