

Which materials are compatible with μ -Slides?

ibidi 社マイクロスライドの薬剤耐性

■耐性がある有機溶剤

Hydrogenperoxide, Sodium hydroxide, Sulfuric acid	+++
DMSO, Acetonitrile	+++
Acetone, Methanol, Ethanol, Isopropanol	+++
Formaldehyde	+++
Ethyl acetate	++

■耐性がない有機溶剤

Pentane, Heptane, Hexane	---
Oleic acid	---
Toluene, Benzene	---
Pyridine / Lutidin	---
THF, Chloroform	---
Modeling clay / Plasticine	---
Mineral oil	---